



⑤ Diagnosis EKG Cables & Accessories


Specification

Cable material	TPU Jacket, Light gray
Connector & Yoke	Nylon plug & Golden plated pins
Trunk Cable length	2m
Lead wire length	1.2m(limb lead);0.8m(Chest lead)

Compatibility

Major brands EKG diagnosis system

Application Diagram





Product name	Photo	REF# Example	Reference to REF# Descriptions
One piece Fixed Cable		E10R-HP-B/I	E: EKG diagnostic cable or lead wire 10: 10-lead(12=12 lead) R: built in resistors for protection(default: no resistors) HP: Philips/HP,factory/model code I: Connector configure; B: banana end(N= needle, P=grabber, S=snap) I: IEC code (Default is AHA code)
Trunk Cable		E10R-MQ-I	E: EKG diagnostic cable or lead wire 10: 10-lead(12=12 lead) R: built in resistors for protection(default: no resistors) MQ: GE Marquette, factory/model code I: Connector configure (default: only one connector available) I: IEC code (Default is AHA code)
Lead Wire		MQ10-L-B-I	MQ: Marquette, factory/model code 10: 10: 10-lead(12=12 lead) L: lead wires B: banana end(N= needle, P=grabber, S=snap) I: IEC code (Default is AHA code)



All one piece cable or lead wires are available with banana end, needle end, grabber end and snap end

Code	AHA(A)	Standard IEC (I)	Alternate IEC (IH)
4mm Banana end			
3mm Needle end			
Grabber End			
Snap End			

5.2 EKG Trunk Cable



OEM & Models	EKG machine Plug	Lead#	REF#	Used with lead wire set/GE Multi-link set	Price code
Burdick					
EK10-E350		10	E10-BK1 E10R-BK1	10-lead, chest /0.9m;Limb/1.2m MQ10-LB,MQ10-LN,MQ10-LP,MQ10-LS or MQB10-LB,...	D3
		10	E10-BK1-I E10R-BK1-I	10-lead, chest /0.9m;Limb/1.2m MQ10-LB-I,MQ10-LN-I,MQ10-LP-I,MQ10-LS-I or MQB10-LB, ... (See the next page of 10-lead wire sets for detail)	D3
Burdick Quinton Eclipse 4,8,400,800, 850,850i,plus, Eclipse LE II and ATRIA PLUS		10	E10-BK2 E10R-BK2	Same as above	D3
		10	E10-BK2-I E10R-BK2-I	Same as above	D3
Fukuda Denshi					
FX101, FX102F, FX328U, FCP2201		10	E10-FD1 E10R-FD1	Same as above	D3
		10	E10-FD1-I E10R-FD1-I	Same as above	D3
Fukuda ME					
KP-500, KP-500D		10	E10-FD2 E10R-FD2	Same as above	D3
		10	E10-FD2-I E10R-FD2-I	Same as above	D3

OEM& Models	EKG machine Plug	Lead#	REF#	Used with lead wire set/GE Multi-link set	Price code
GE-Marquette					
MAC 500, MAC 1100, MAC 1200		10	E10-MQ E10R-MQ	Same as above	D3
		10	E10-MQ-I E10R-MQ-I	Same as above	D3
ELI-200, ELI-50		10	E10-MT E10R-MT	Same as above	D3
		10	E10-MT-I E10R-MT-I	Same as above	D3
Hellige / Siemens Hormann / Bosch					
		10	E10-HL E10R-HL	Same as above	D3
		10	E10-HL-I E10R-HL-I	Same as above	D3
Kenz					
103,106 CARDIOLINE (Italy)		10	E10-KZ1 E10R-KZ1	Same as above	D3
		10	E10-KZ1-I E10R-KZ1-I	Same as above	D3
ECG 108/110 Cardico 302/601/ 1210/1211		10	E10R-KZ2	Same as above	D3
		10	E10R-KZ2-I	Same as above	D3
Nihon Kohden					
CaduoFax 6151, 6851,6353 CaduoFax 8110		10	E10-NK1 E10R-NK1	Same as above	D3
		10	E10-NK1-I E10R-NK1-I	Same as above	D3
CaduoFax Q ECG-9110K ECG-9130K Doing Juang Cardioline (New)		10	E10-NK2 E10R-NK2	Same as above	D3
		10	E10-NK2-I E10R-NK2-I	Same as above	D3
Philips/HP					
M1770A,M1771A, M1772A Pagewrite100/200/ 200i		10	E10-HP E10R-HP	Same as above	D3
		10	E10-HP-I E10R-HP-I	Same as above	D3
Schiller					
AT1, AT2, AT2 Plus,AT4 Resistor 10K Long Screw		10	E10-SH1 E10R-SH1	Same as above	D3
		10	E10-SH1-I E10R-SH1-I	Same as above	D3
AT3,AT6,CS6,AT5, AT10,AT60 Short Screw		10	E10-SH2 E10R-SH2	Same as above	D3
		10	E10-SH2-I E10R-SH2-I	Same as above	D3
GE Healthcare-Resting sys.					
CAM14 coiled patient Cable Used with GE MAC 5000 resting system with CAM14 Acquisition Module		10	ECAM-GE14-S 2016560-001 1.4m+0.3m ECAM-GE14 2016560-003 4.5m+0.65m	Used with CAM 14 patient lead wires on next page.	D8
		10	EAM-GE4-S 700044-203 1.3m+0.25m EAM-GE4 700044-202 4.6m+0.25m	Used with CAM 4 patient lead wires on next page.	D8

5.3 EKG Lead Wire


Pictures	OEM & Models	Electrode end	REF# (AHA)	REF# (IEC)	Pieces/Set	Price code
GE Multi-Link 10-Lead-integrated leads						
	Banana End 900112-001(AHA) 900112-002(IEC)	4mm banana	MQ10-LB 900112-001	MQ10-LB-I 900112-002	10	D1
		3mm needle	MQ10-LN	MQ10-LN-I	10	D1
		Grabber	MQ10-LP	MQ10-LP-I	10	D1
		Snap	MQ10-LS	MQ10-LS-I	10	D1
GE CAM 14 patient lead wire						
		4mm banana	MQB10-LB	MQB10-LB-I	10	D1
		3mm needle	MQB10-LN	MQB10-LN-I	10	D1
		Grabber	MQB10-LP	MQB10-LP-I	10	D1
		Snap	MQB10-LS	MQB10-LS-I	10	D1
	420101-002	CAM/Universal	MQB10-AD	MQB10-AD-I	10	D1
	420101-001 Used with 900178-001/ 900178-002 or Unimed Adapter as electrodes end (see below)	CAM/Universal	MQB14-AD	MQB14-AD-I	14	D3
	2001925-005(66cm)	CAM/Universal	X-MQB-66DF	66cm, no color coded	10	C8
	2001925-002(74cm)	CAM/Universal	X-MQB-80DF	80cm, no color coded	10	C8
	2001925-006(90cm)	CAM/Universal	X-MQB-90DF	90cm, no color coded	10	C8
	2001925-003(102cm)	CAM/Universal	X-MQB-100DF	100cm, no color coded	10	C8
	2001925-004(130cm)	CAM/Universal	X-MQB-120DF	120cm, no color coded	10	C8
	2001925-003(102cm)	CAM/Universal	X-MQB-80DF	80cm, no color coded	10	C8
2001925-004(130cm)	CAM/Universal	X-MQB-120DF	120cm, no color coded	10	C8	
GE AM4/AM5 patient lead wires						
	Banana End 900112-001(AHA) 900112-002(IEC)	4mm banana	MQD10-LB	MQD10-LB-I	10	D1
		3mm needle	MQD10-LN	MQD10-LN-I	10	D1
		Grabber	MQD10-LP	MQD10-LP-I	10	D1
		Snap	MQD10-LS	MQD10-LS-I	10	D1
	900177-001 Used with 900178-001/ 900178-002 or Unimed Adapter as electrodes end (see below)	CAM/Universal	MQD12-AD	MQD12-AD-I	12	D3
	406581-004(90cm)	CAM/Universal	X-MQD-80DF	80cm, no color coded	10	C8
	406581-006(130cm)	CAM/Universal	X-MQD-120DF	120cm, no color coded	10	C8
GE CAM14/AM4/AM5 Electrode Adapter						
	900178-001,AHA 900178-101,IEC 900179-201,AHA 900178-201,IEC For the above leads set or GE original CAM14/AM4/AM5	4mm banana, 10 pcs/set	AD10-DB	AD10-DB-I	10	B8
		4mm banana, 14 pcs/set	AD14-DB	AD14-DB-I	14	C3
		4mm banana, 12 pcs/set	AD12-DB	AD12-DB-I	12	C2

Pictures	OEM & Models	Electrode end	REF# (AHA)	REF# (IEC)	Pieces/Set	Price code
Philips Lead wire						
	PageWriter Touch: 860284, PageWriter Trim Cardiographs, 860286, 860287, 860288, 860290, 860297	4mm banana	PW10-LB 989803129161	PW10-LB-I 989803129191	10	D1
		3mm needle	PW10-LN	PW10-LN-I	10	D1
		Grabber	PW10-LP	PW10-LP-I	10	D1
		Snap	PW10-LS	PW10-LS-I	10	D1
	989803151631/989803151641	4mm banana	PWC10-LB	PWC10-LB-I	10	D1
		3mm needle	PWC10-LN	PWC10-LN-I	10	D1
		Grabber	PWC10-LP	PWC10-LP-I	10	D1
		Snap	PWC10-LS	PWC10-LS-I	10	D1
QUINTON Lead wire						
	Q-710 Stress Test System 031257-002	4mm banana	QT10-LB	QT10-LB-I	10	D1
		3mm needle	QT10-LN	QT10-LN-I	10	D1
		Grabber	QT10-LP	QT10-LP-I	10	D1
		Snap	QT10-LS	QT10-LS-I	10	D1
Mortara Lead wire						
	9293-046-60 ELI150C ELI230 ELI250C ELI350	4mm banana	MT10-LB	MT10-LB-I	10	D1
		3mm needle	MT10-LN	MT10-LN-I	10	D1
		Grabber	MT10-LP	MT10-LP-I	10	D1
		Snap	MT10-LS	MT10-LS-I	10	D1
Norav (AA style) Lead wire						
		4mm banana	AA10-LB	AA10-LB-I	10	D1
		3mm needle	AA10-LN	AA10-LN-I	10	D1
		Grabber	AA10-LP	AA10-LP-I	10	D1
		Snap	AA10-LS	AA10-LS-I	10	D1
Welch Allyn Lead wire						
	CP100/CP200 401129	4mm banana	WA10-LB	WA10-LB-I	10	D1
		3mm needle	WA10-LN	WA10-LN-I	10	D1
		Grabber	WA10-LP	WA10-LP-I	10	D1
		Snap	WA10-LS	WA10-LS-I	10	D1

5.4 EKG Accessories 附件



ECG Electrodes 心电电极	REF# 型号	Package 包装	Price code 价格代码	REF to Figures 图例
Adult Limb Clamps/IEC, universal 成人四肢夹/欧标, 通用	UU-LC-A	4 pcs/box	B6	A
Adult Limb Clamps/AHA, all clamps in green 成人四肢夹/美标, 绿色	UU-LC-AA	4 pcs/box	B6	E
Pediatric Limb Clamps/IEC, universal 儿童四肢夹/欧标, 通用	UU-LC-P	4 pcs/box	B6	B
Pediatric Limb Clamps/AHA, all clamps in green 儿童四肢夹/美标, 绿色	UU-LC-PA	4 pcs/box	B6	/
Adult up Limb Clamps, for snap electrode, 3-lead 成人上肢夹, 接扣式电极, 3极	UU-DLC-A	2 pcs/box	B6	F
Adult suction electrodes, Universal 成人胸吸球, 通用	UU-SE-A	6 pcs/box	B8	C
Pediatric suction electrodes, Universal 儿童胸吸球, 通用	UU-SE-P	6 pcs/box	B8	D
Adult blue bulb 成人蓝色吸球	Y05-F027	6 pcs/box	A8	G
Pediatric green bulb 儿童绿色吸球	Y05-F022	6 pcs/box	A8	H



ECG Electrodes 心电电极	REF# 型号	Package 包装	Price code 价格代码	REF to Figures 图例
* Flat clip/red Adapter/w 4mm banana 平式转换头/红色/4mm 香蕉头	AD-BF-10	10 pcs/bag	B7	A
* Flat clip/blue adapter/w 3mm needle 平式转换头/蓝色/3mm 针头	AD-NF-10	10 pcs/bag	B7	B
Alligator clip/red Adapter/w 4mm banana 鳄鱼夹/红色/4mm 香蕉头	AD-BA-10	10 pcs/bag	B7	C
Alligator clip/blue adapter/w 3mm needle 鳄鱼夹/蓝色/3mm 针头	AD-NA-10	10 pcs/bag	B7	D
Adapter/4mm with grabber end 转接头/4mm 钳式	AD-BP-10	10 pcs/bag	B8	E
Adapter/4mm with snap end 转接头/4mm 扣式	AD-BS-10	10 pcs/bag	B8	F
Wide mouth Multi-function adapter/4mm with banana end 多功能心电转换钳/4mm 香蕉头	AD-BW-10	10 pcs/bag	B8	I
Wide mouth Multi-function adapter/3mm with banana end 多功能心电转换钳/3mm 香蕉头	AD-NW-10	10 pcs/bag	B8	J
Adapter wire/4mm with grabber end, 0.3m 转接线/4mm/香蕉头>>钳式	AD-B-P	10 pcs/bag	B8	G
Adapter wire/3mm with grabber end, 0.3m 转接线/3mm/针式>>钳式	AD-N-P	10 pcs/bag	B8	G
Adapter wire/4mm with snap end, 0.3m 转接线/4mm/香蕉头>扣式	AD-B-S	10 pcs/bag	B8	H
Adapter wire/3mm with snap end, 0.3m 转接线/3mm/针式>扣式	AD-N-S	10 pcs/bag	B8	H

Items indicated "★" are not sold in the USA.
带*的产品均不在美国区域销售