

4 Accessories



REF#	Product Name	Purpose of Use	Package	Price	Figures
Y04-E041	Ear sensor hanger	Used with Unimed ear clip sensor U900 series or Other brands ear sensor(Gray)	1 pcs	A5	A
UEC-01	Adult ear clip Or as animal small size lingua clip	Used with D-YS or Unimed U800 sensor(Black)	1 pcs	A5	B
UVC-01	Animal lingua clip,Big size	Used with D-YS or Unimed U800 sensor(Gray)	1 pcs	A8	C
UVC-02	Animal lingua clip,Middle size	Used with D-YS or Unimed U800 sensor(Gray)	1 pcs	A7	D
Y04-E019	Reusable Silicon Y sensor wrap tape	Used with Nellcor D-YS or Unimed U800 series sensors	1 pcs	A6	E
E023-12	Reusable 3 holes foam tape For All patient sizes	Used with Nellcor D-YS or Unimed U800 /U600 series sensors	12 pcs Per bag	B5/bag	F
E024-25	Disposable 3M adhesive textile tape For Adult and neonate	Used with Nellcor D-YS or Unimed U800 series sensors	25 pcs Per bag	B7/bag	G
E025-25	Disposable 3M adhesive textile tape for Pediatric and infant	Used with Nellcor D-YS or Unimed U800 series sensors	25 pcs Per bag	B7/bag	H
T002-12	Velcro foam wrap tape	To fix the sensor cable onto the patient wrist	12 pcs per bag	C1.4/bag	I