



1 Reusable SpO₂ Sensor

Specification

Sensor Connector	OEM compatible, Nylon adapter with gold plated pins
Cable material & color	TPU jacket; light gray or blue
Cable length & diameter	1.1m/1.5m/3.0m; ϕ3.20mm/4.00mm or 5.00x2.50mmtwo line cable
Sensor house	ABS shell(finger clip); Silicon rubber boots(soft and multi-site)

Compatible With Major Patient Monitor Brands

- Nellcor

Ohmeda

Philips

Marquette

Spacelabs

Novamatrix
- Criticare

BCI


Datascope

Datex

Nihon Kohden

Palco

.....



Series系列	U400	U100	U400S	U100S	U200S	U300
REF to Fig	A	B	C	D	E	F
Name	Adult Finger	Pediatric Finger	Adult soft tip	Pediatric Soft tip	Infant Soft tip	Neonate Wrap
Patient size	+40kg	15-40kg	+40kg	15-40kg	3-15kg	1-3kg
Cable color	Gray Cable	Gray Cable	Blue Cable	Gray two line	Gray two line	Gray two line
Series系列	U800	U900	U600	U600B	U800E	U800V
REF to Fig	G	H	I	J	K	L
Name	Multi-site	Adult Ear clip	Adult/Neo. Wrap	Pedia/Inf. Wrap	Multi-site	Forehead/Rectal
Patient size	All patients	+40kg	+40kg or -3kg	15-40kg	All patients	Adult
Cable color	Gray two line		Gray two line	Gray two line	Sensormate Used with DPT tapes	Gray Cable

How to Identify the REF# of Unimed SpO₂ Sensors



Reference Number Rule Description

U403S-01XE	
U	Unimed brand
4	adult sensor (1: pediatric; 2: infant;3: neonate; 5: disposable; 6 wrap; 7: ext-cable; 8: multi-site;9: ear clip)
03	1.1m cable length (03: 1.1m;05:1.5m;08: 2.2m;10: 3.0m)
S	Soft sensor (B:U style shell; Default: shell sensor or neonate wrap, alternative wrap, multi-site)
01	Nellcor code (02:Ohmeda;03:Novametrix; 04: Spacelabs; 05: Criticare; 06: BCI...)
X	Oximax technology(M: M technology, D: digital technology; Default: standard
E	Sensormate sensor (Default: non-Sensormate sensor)

Cross Reference

Ending "X" can be replaced by "XP" or "P" for fully compatible with all oximax oximetry

OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.
3F Medical (China)								
IRIS 300/IRIS 400 Phoebe PMS8300/8310 PhoebeProPMS8300/ 8310PMS8110A			U410-189	Adult finger /3.0m	C7	No need		
			U410S-189	Adult soft /3.0m	C8			
			U810-189	Multi-site Y/3.0m	C7			
Adecon (China)								
MCF-21L, DK-8000D,DK-8000L			U410-120	Adult finger /3.0m	C7	No need		
			U410S-120	Adult soft /3.0m	C8			
			U810-120	Multi-site Y/3.0m	C7			
Anke (China)								
ASC-553			U410-37	Adult finger /3.0m	C7	No need		
			U410S-37	Adult soft成/3.0m	C8			
			U810-37	Multi-site Y/3.0m	C7			
ASC-553 (*Oximax tech)			U410-78P	Adult finger /3.0m	C8	No need		
			U410S-78P	Adult soft /3.0m	D1			
			U810-78P	Multi-site Y/3.0m	C8			
ASC-553A			U410-82	Adult finger /3.0m	C7	No need		
			U410S-82	Adult soft /3.0m	C8			
			U810-82	Multi-site Y/3.0m	C7			
ASC-553A (*Ohmeda Tech)			U410-107	Adult finger /3.0m	C7	No need		
			U410S-107	Adult soft /3.0m	C8			
			U810-107	Multi-site Y/3.0m	C7			

OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.
Artema/S&W								
Diascope,8000, Athena (*Non-Oximax)		DS100A	U403-01	Adult finger/1.1m	C3	U708-80		
		DS100A	U403S-01	Adult soft /1.1m	C5			
		D-YS	U803-01	Multi-site Y/1.1m	C3			
			U410-80	Adult finger /3.0m	C8	No need		
			U410S-80	Adult soft /3.0m	D1			
			U810-80	Multi-site Y/3.0m	C8			
S&W			U410-19	Adult finger /3.0m	C8	No need		
			U410S-19	Adult soft成 /3.0m	D1			
			U810-19	Multi-site Y/3.0m	C8			
Baxter/Simed								
ASAT, S-50/100/100a			U410-12	Adult finger /3.0m	C8	No need		
			U410S-12	Adult soft /3.0m	D1			
			U810-12	Multi-site Y/3.0m	C8			
Monarch/Explorer (*Oximax)		DS100A	U403-01	Adult finger/1.1m	C3	U708-01		
		DS100A	U403S-01	Adult soft/1.1m	C5			
		D-YS	U803-01	Multi-site Y/1.1m	C3			
BCI								
3300,3301,3302, 3303,3304,3401 (*BCI)		3044	U403-06	Adult finge /1.1m	C3	Depends U708-06		
			U103-06	Pediatric finger /1.1m	C5			
			U403S-06	Adult soft /1.1m	C3			
			U103S-06	Pediatric soft /1.1m	C5			
			U203S-06	Infant soft /1.1m	C5			
			U303-06	Neonate wrap /1.1m	C5			
		3043	U803-06	Multi-site Y/1.1m	C3			
			U603-06	Adult/Neo Wrap/1.1m	C3			
			U603B-06	Pedia/infant wrap/1.1m	C3			
		3078	U903-06	Adult ear clip /1.1m	C5			
3300,3301 6100,9100 Advisor,Autocorr, MiniToor, Mini Torr plus (BCI*)		See above *	See the above sensor*	1.1m	See above sensors*	U708-76		
			U410-76	Adult finger /3.0m	C7	No need		
			U410S-76	Adult soft /3.0m	C8			
			U810-76	Multi-site Y/3.0m	C7			

OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.
Biosys								
BPM200 BPM700 Guardian		DS100A	U403-01	Adult finger /1.1m	C3	U708-01		
		DS100A	U403S-01	Adult soft /1.1m	C5			
		D-YS	U803-01	Multi-site Y/1.1m	C3			
			U410-27	Adult finger /3.0m	C8	No need		
			U410S-27	Adult soft /3.0m	D1			
			U810-27	Multi-site Y/3.0m	C8			
Biolight(China)								
M6,M12,M69			U410-47	Adult finger /3.0m	C7	No need		
			U410S-47	Adult soft/3.0m	C8			
			U810-47	Multi-site Y/3.0m	C7			
M66,M69			U410-57	Adult finger /3.0m	C8	No need		
			U410S-57	Adult soft /3.0m	D1			
			U810-57	Multi-site Y/3.0m	C8			
(*Oximax)			U410-51P	Adult finger /3.0m	C8	No need		
			U410S-51P	Adult soft/3.0m	D1			
			U810-51P	Multi-site Y/3.0m	C8			
M700 (*Digital Tech)			U403-62D	Adult finger /1.1m	C8	No need		
			U403S-62D	Adult soft /1.1m	D1			
			U803-62D	Multi-site Y/1.1m	C8			
M8000,M9000,M8000A M9000A,BTD-352A, new M69 (*Digital Tech)			U410-62D	Adult finger /3.0m	D1	No need		
			U410S-62D	Adult soft /3.0m	D2			
			U810-62D	Multi-site Y/3.0m	D1			
M900,M7000,M9500 A8 Modular (*Digital Tech)			U410-68D	Adult finger /3.0m	D1	No need		
			U410S-68D	Adult soft /3.0m	D2			
			U810-68D	Multi-site Y/3.0m	D1			
Anyview A2E Anyview A3 V6vet			U410-205D	Adult finger /3.0m	C8	No need		
			U410S-205D	Adult soft/3.0m	D1			
			U810-205D	Multi-site Y/3.0m	C8			
M800 handled vet			U403-62D	Adult finger /1.1m	C8	U708-205		
			U403S-62D	Adult soft/1.1m	D1			
			U803-62D	Multi-site Y/1.1m	C8			
Anyview A3 (*Masimo tech)			U410-87	Adult finger /3.0m	D2	No need		
			U410S-87	Adult soft/3.0m	D3			
			U810-87	Multi-site Y/3.0m	D2			
Bionet								
		DS100A	U403-01	Adult finger /1.1m	C3	U708-58		
		DS100A	U403S-01	Adult soft /1.1m	C5			
		D-YS	U803-01	Multi-site Y/1.1m	C3			
BM3/BM5			U410-58	Adult finger /3.0m	C7	No need		
			U410S-58	Adult soft /3.0m	C8			
			U810-58	Multi-site Y/3.0m	C7			
BM3/BM5(New)			U410-108	Adult finger /3.0m	C7	No need		
			U410S-108	Adult soft /3.0m	C8			
			U810-108	Multi-site Y/3.0m	C7			

OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.
Biocare (China)								
BP-12B,BP-12C (*Digital Tech)			U410-61D	Adult finger /3.0m	C8	No need		
			U410S-61D	Adult soft/3.0m	D1			
			U810-61D	Multi-site Y/3.0m	C8			
BM-9000,BM-9000S FM-801,FM-801V (*Digital Tech)			U410-62D	Adult finger /3.0m	D1	No need		
			U410S-62D	Adult soft/3.0m	D2			
			U810-62D	Multi-site Y/3.0m	D1			
BM-9000			U410-29	Adult finger /3.0m	C7	No need		
			U410S-29	Adult soft/3.0m	C8			
			U810-29	Multi-site Y/3.0m	C7			
IM15			U410-223	Adult finger /3.0m	C8	No need		
			U410S-223	Adult soft/3.0m	D1			
			U810-223	Multi-site Y/3.0m	C8			
PM-900			U410-143	Adult finger /3.0m	C7	No need		
			U410S-143	Adult soft/3.0m	C8			
			U810-143	Multi-site Y/3.0m	C7			
CAS Medical								
Oscillomate 9002 9303,9302S (*Nonin)		8000AA	U403-08	Adult finger /1.1m	C3	U708-01		
			U403S-08	Adult soft/1.1m	C5			
			U803-08	Multi-site Y/1.1m	C3			
740-750(M) (*M tech)		1863/ LNCS DCI	U403-49R	Adult finger /1.1m	C4	U708M-15R		
		1863/ LNCS DCI	U403S-49R	Adult soft/1.1m	C6			
		2258/ LNCS YI	U803-49R	Multi-site Y/1.1m	C4			
		LNOP DC-12	U410-15R	Adult finger /3.0m	C8	No need		
		LNOP DC-12	U410S-15R	Adult soft/3.0m	D1			
			U810-15R	Multi-site Y/3.0m	C8			
Charmcare								
CX100 (*Non-Oximax)		DS100A	U403-01	Adult finger /1.1m	C3	U708-01		
		DS100A	U403S-01	Adult soft/1.1m	C5			
		D-YS	U803-01	Multi-site Y/1.1m	C3			
ChoiceMed超思(China)								
MMED3001		DS100A	U403-01	Adult finger /1.1m	C3			
		DS100A	U403S-01	Adult soft/1.1m	C5			
		D-YS	U803-01	Multi-site Y/1.1m	C3			
MMED100H, MMED300A, MMED6000DP (*Digital tech)			U410-84D	Adult finger /3.0m	D1	No need		
			U410S-84D	Adult soft/3.0m	D2			
			U810-84D	Multi-site Y/3.0m	D1			
MMED6000DP (*BCI Tech)			U410-96	Adult finger /3.0m	C7	No need		
			U410S-96	Adult soft/3.0m	C8			
			U810-96	Multi-site Y/3.0m	C7			
MMED6000DP			U410-37	Adult finger /3.0m	C7	No need		
			U410S-37	Adult soft/3.0m	C8			
			U810-37	Multi-site Y/3.0m	C7			
ChoiceMed(China)								

OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.			
ChoiceMed(China)											
MMED3001,MD300K MA300K2,MD300M, MMED6000DP-M3 (*BCI Tech)			U403-06	Adult finger /1.1m	C3	No need					
			U403S-06	Adult soft /1.1m	C5						
			U803-06	Multi-site Y/1.1m	C3						
Contec康泰(China)											
			U410-29	Adult finger /3.0m	C7	No need					
			U410S-29	Adult soft /3.0m	C8						
			U810-29	Multi-site Y/3.0m	C7						
CMS60C/CW/D/DW,CMS60, PM50,PM60D,CSM5000, CMS5000C,CMS5100, CMS6000,CMS7000/ 8000/9000(*Digital Tech)			U410-46D	Adult finger /3.0m	C8	No need					
			U410S-46D	Adult soft /3.0m	D2						
			U810-46D	Multi-site Y/3.0m	C8						
CMS50F/CMS50FW			U401-149	Adult finger /0.3m	C7	No need					
			U401S-149	Adult soft /0.3m	C8						
Colin/Omron											
BP306,BP88,BP508 (Nellcor Non-Oximax) M570		DS100A	U403-01	Adult finger /1.1m	C3	U708-01					
		DS100A	U403S-01	Adult soft /1.1m	C5						
		D-YS	U803-01	Multi-site Y/1.1m	C3						
BP88S (*Non-Oximax)		See the above	See the above	See the above	See the above	U708-69					
			U410-69	Adult finger /3.0m	C8				No need		
			U410S-69	Adult soft /3.0m	D1						
			U810-69	Multi-site Y/3.0m	C8						
BX 10NE (* Oximax)		DS100A	U403-01XP	Adult finger /1.1m	C5	U710-70					
		DS100A	U403S-01XP	Adult soft /1.1m	C7						
		D-YS	U803-01XP	Multi-site Y/1.1m	C5						
BX 10MA (* M tech)		1863/ LNCS DCI	U403-49R	Adult finger /1.1m	C4	U708M-15R					
		1863/ LNCS DCI	U403S-49R	Adult soft /1.1m	C6						
		2258/ LNCS YI	U803-49R	Multi-site Y/1.1m	C4						
		LNOP DC-12	U410-15R	Adult finger /3.0m	C8	No need					
		LNOP DC-12	U410S-15R	Adult soft /3.0m	D1						
			U810-15R	Multi-site Y/3.0m	C8						
Corpuls											
08/16 (*BCI)		3044	U403-06	Adult finger /1.1m	C3	NA					
			U403S-06	Adult soft /1.1m	C5						
		3043	U803-06	Multi-site Y/1.1m	C3						
08/16(M) (* Non-Oximax)		DS100A	U403-01	Adult finger /1.1m	C3	NA					
		DS100A	U403S-01	Adult soft /1.1m	C5						
		D-YS	U803-01	Multi-site Y/1.1m	C3						
08/16(M) M Tech			U410-152M	Adult finger /3.0m	C8	No need					
			U410S-152M	Adult soft /3.0m	D1						
			U810-152M	Multi-site Y/3.0m	C8						

OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.
Comen (China)								
Start-8000A/8000B/8000C Star-8000D,8000H Star-5000B/5000C,C50 (*Digital Tech)			U410-66D	Adult finger /3.0m	C8	No need		
			U410S-66D	Adult soft /3.0m	D1			
			U810-66D	Multi-site Y/3.0m	C8			
Start-8000E			U410-183	Adult finger /3.0m	C7	No need		
			U410S-183	Adult soft /3.0m	C8			
			U810-183	Multi-site Y/3.0m	C7			
C60,C80 (*Oximax Tech)			U410-187P	Adult finger /3.0m	D1	No need		
			U410S-187P	Adult soft /3.0m	D2			
			U810-187P	Multi-site Y/3.0m	D1			
C100/C100A/C21/C30 C50/C70/C80/NC3 NC5/NC8/NC12/NC19			U410-208D	Adult finger /3.0m	C8	No need		
			U410S-208D	Adult soft /3.0m	D1			
			U810-208D	Multi-site Y/3.0m	C8			
Old C60 Masimo Rainbow			U410-202M	Adult finger /3.0m	D2	No need		
			U410S-202M	Adult soft /3.0m	D3			
			U810-202M	Multi-site Y/3.0m	D2			
C20,C60,C80,C90 (*Nellcor OxiMax Tech)			U410-208P	Adult finger /3.0m	D1	No need		
			U410S-208P	Adult soft /3.0m	D2			
			U810-208P	Multi-site Y/3.0m	D1			
C20,C60,C80,C90 (Masimo Rainbow)			U410-208MR	Adult finger /3.0m	D2	No need		
			U410S-208MR	Adult soft /3.0m	D3			
			U810-208MR	Multi-site Y/3.0m	D2			
C60,C80 (*Masimo Tech)			U410-187M	Adult finger /3.0m	D2	No need		
			U410S-187M	Adult soft /3.0m	D3			
			U810-187M	Multi-site Y/3.0m	D2			
Creative(China)								
UP-6000,PC-9000 (*Non-Oximax)		DS100A	U403-01	Adult finger /1.1m	C3	U708-29		
		DS100A	U403S-01	Adult soft /1.1m	C5			
		D-YS	U803-01	Multi-site Y/1.1m	C3			
			U410-29	Adult finger /3.0m	C7	No need		
			U410S-29	Adult soft /3.0m	C8			
UP-7000,UP-9000,PC-900, PG-900A,PG-900S PG-900SNPC-3000, PG-5000A/B,PC-9000B, PG-66A, PC-66B, PG-66V(*Digital tech)			U410-61D	Adult finger /3.0m	C8	No need		
			U410S-61D	Adult soft /3.0m	D1			
			U810-61D	Multi-site Y/3.0m	C8			

OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.
Criticare (CSI)								
510,503,504, 507,508,581, 1100,8100,Poet, Scholar (CSI*)		934-10LN	U410-05	Adult finger /3.0m	C8	No need		
			U110-05	Pediatric finger /3.0m	D1			
		975AD-10	U410S-05	Adult soft /3.0m	D1			
			U110S-05	Pediatric soft /3.0m	D1			
			U210S-05	Infant soft /3.0m	D1			
			U310-05	Neonate wrap /3.0m	D1			
			U610-05	Adult/Neo wrap /3.0m	C8			
			U610B-05	Pedia/infant wrap/3.0m	C8			
			U810-05	Multi-site Y/3.0m	C8			
			U910-05	Adult ear clip /3.0m	D1			
			U403-05	Adult finger /1.1m	C3	U708-05		
			U103-05	Pediatric finger /1.1m	C4			
			U403S-05	Adult soft /1.1m	C4			
			U103S-05	Pediatric soft /1.1m	C4			
			U203S-05	Infant soft /1.1m	C4			
			U303-75	Neonate wrap /1.1m	C4			
			U603-05	Multi-site Y/1.1m	C3			
			U603B-05	Adult/Neo/1.1m	C3			
			U803-05	Pedia/infant wrap/1.1m	C4			
			U903-05	Adult ear clip /1.1m	C4			
503SPOT,504Plus, 504DX,506DX (*CSI)		934-10DN	U410-75	Adult finger /3.0m	C7	No need		
			U110-75	Pediatric finger /3.0m	C8			
			U410S-75	Adult soft /3.0m	C8			
			U110S-75	Pediatric soft /3.0m	C8			
			U210S-75	Infant soft /3.0m	C8			
			U310-75	Neonate wrap /3.0m	C8			
			U610-75	Adult/Neo wrap/3.0m	C7			
			U610B-75	Pedia/infant wrap/3.0m	C7			
			U810-75	Multi-site Y/3.0m	C7			
			U910-75	Adult ear clip /3.0m	C8			
503SPOT,503DX,504Plus, 504DX,506DX,506DXNP (*CSI)			U410-75N	Adult finger/3.0m	C7	No need		
			U410S-75N	Adult soft/3.0m	C8			
			U810-75N	Multi-site Y/3.0m	C7			
Digicare								
Digimax5500 (*BCI)		3044	U403-06	Adult finger /1.1m	C3	U708-76		
			U403S-06	Adult soft /1.1m	C5			
		3043	U803-06	Multi-site Y/1.1m	C3	No need		
			U410-76	Adult finger /3.0m	C7			
			U410S-76	Adult soft /3.0m	C8			
			U810-76	Multi-site Y/3.0m	C7	No need		
			U410-52	Adult finger /3.0m	C7			
			U410S-52	Adult soft /3.0m	C8			
			U810-52	Multi-site Y/3.0m	C7	No need		


OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.
Datascope								
Accusat, Accutor 3/4 sat Passport, Accutor plus (*Datascope)			U403-07	Adult finger /1.1m	C4	U708-07		
			U103-07	Pediatric finger /1.1m	C6			
			U403S-07	Adult soft /1.1m	C6			
			U103S-07	Pediatric soft /1.1m	C6			
			U203S-07	Infant soft /1.1m	C6			
			U303-07	Neonate wrap /1.1m	C6			
			U603-07	Adult/Neo wrap/1.1m	C5			
			U603B-07	Pedia/infant wrap/1.1m	C5			
			U803-07	Multi-site Y/1.1m	C5			
		U903-07	Adult ear clip /1.1m	C4				
			U410-07	Adult finger /3.0m	C8	No need		
			U410S-07	Adult soft /3.0m	D1			
			U810-07	Multi-site Y/3.0m	C8			
Accutor plus (*M Tech)		1863/ LNCS DCI	U403-49R	Adult finger /1.1m	C4	U708M-15R		
		1863/ LNCS DCI	U403S-49R	Adult soft /1.1m	C6			
		2258/ LNCS YI	U803-49R	Multi-site Y/1.1m	C4			
		LNOP DC-12	U410-15	Adult finger /3.0m	C8	No need		
		LNOP DC-12	U410S-15	Adult soft /3.0m	D1			
			U810-15	Multi-site Y/3.0m	C8			
PassportII, Expert Pssport XG, Accutor plus (*Non-Oximax)		DS100A	U403-01	Adult finger /1.1m	C3	U708-69		
		DS100A	U403S-01	Adult soft /1.1m	C5			
		D-YS	U803-01	Multi-site Y/1.1m	C3			
Accutor-Trio (*Non-Oximax)		See the above	See the above	See the above	See the above	U708-01		
Trio Nouvelle Generation (* Oximax)		DS100A	U403-01XP	Adult finger /1.1m	C5	U710-70		
		DS100A	U403S-01XP	Adult soft /1.1m	C7			
		D-YS	U803-01XP	Multi-site Y/1.1m	C5			
Dima (Italia)								
(*Non-Oximax)		DS100A	U403-01	Adult finger /1.1m	C3	U708-01		
		DS100A	U403S-01	Adult soft /1.1m	C5			
		D-YS	U803-01	Multi-site Y/1.1m	C3			
Dixon								
D6 D8			U410-205	Adult finger /3.0m	C8	No need		
		U410S-205	Adult soft/3.0m	D1				
		U810-205	Multi-site Y/3.0m	C8				
			U403-62D	Adult finger /1.1m	C8	U708-208		
			U403S-62D	Adult soft/1.1m	D1			
			U803-62D	Multi-site Y/1.1m	C8			

OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext -Cable	Cable Plug	Sensor Conn.	
Dixon									
D2,D3,D5,D6,D8 DS3,DS5,DS7 (Digital)			U410-205D	Adult finger /3.0m	C8	No need			
			U410S-205D	Adult soft/3.0m	D1				
			U810-205D	Multi-site Y/3.0m	C8				
M800 handled vet			U403-62D	Adult finger /1.1m	C8	No need			
			U403S-62D	Adult soft/1.1m	D1				
			U803-62D	Multi-site Y/1.1m	C8				
DS3,DS5,DS7 (Masimo Rainbow)			U410-87	Adult finger /3.0m	D2	No need			
			U410S-87	Adult soft/3.0m	D3				
			U810-87	Multi-site Y/3.0m	D2				
Storm 5000			U403-01	Adult finger /1.1m	C3	No need			
			U403S-01	Adult soft/1.1m	C5				
			U803-01	Multi-site Y/1.1m	C3				
Storm 5300,Storm 5500 Storm 5800,Storm 5900			U410-29	Adult finger /3.0m	C7	No need			
			U410S-29	Adult soft/3.0m	C8				
			U810-29	Multi-site Y/3.0m	C7				
Dixtal(Brazil)-Philips									
Full compatible with Novamatrix Sensas		DX-1NC01-0	U410-113	Adult finger /3.0m	C7	No need			
			U410S-113	Adult soft /3.0m	C8				
			U810-113	Multi-site Y/3.0m	C7				
Dolphin									
(*Digital Tech)			U403-50D	Adult finger /1.1m	C8	U708-50			
			U403S-50D	Adult soft /1.1m	D1				
			U803-50D	Multi-site Y/1.1m	C8				
(*Digital Tech)			U410-50D	Adult finger /3.0m	D1	No need			
			U410S-50D	Adult soft /3.0m	D2				
			U810-50D	Multi-site Y/3.0m	D1				
Drager/Siemens									
Oxipac 2000 (*BCI)		3044	U403-06	Adult finger /1.1m	C3	No need			
			U403S-06	Adult soft /1.1m	C5				
			3043	U803-06	Multi-site Y/1.1m				C3
Vitalert 1000,icero (*Non-Oximax)		DS100A	U403-01	Adult finger /1.1m	C3	U708-28			
			DS100A	U403S-01	Adult soft /1.1m				C5
			D-YS	U803-01	Multi-site Y/1.1m				C3
Narkomed-Oxysat 2 (*Non-Oximax)		See the above	See the above	See the above	See the above	U710-85			
Sirecust, Sirem Spo2 (*Non-Oximax)		See the above	See the above	See the above	See the above	U708-33			
SC6000, 7000,9000 (*Non-Oximax)		See the above	See the above	See the above	See the above	U704-23 U708-23			
Drager/Siemens									

OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.
Drager/Siemens								
Oxytrend (*M tech))		1863/ LNCS DCI	U403-49R	Adult finger /1.1m	C4	U708M-15R		
		1863/ LNCS DCI	U403S-49R	Adult soft /1.1m	C6			
		2258/ LNCS YI	U803-49R	Multi-site Y/1.1m	C4			
		LNOP DC-12	U410-15R	Adult finger /3.0m	C8	No need		
		LNOP DC-12	U410S-15R	Adult soft /3.0m	D1			
			U810-15R	Multi-site Y/3.0m	C8			
(*Oximax Tech)			U410-23P	Adult finger /3.0m	D1	No need		
			U410S-23P	Adult soft/3.0m	D2			
			U810-23P	Multi-site Y/3.0m	D1			
(*M Tech)			U410-23MR	Adult finger /3.0m	C8	No need		
			U410S-23MR	Adult soft/3.0m	D1			
			U810-23MRT	Multi-site Y/3.0m	C8			
Sirecust 1281 Monitor SC 6002XL,SC7000, SC9000XL, (Masimo Rainbow)			U410-23MR	Adult finger /3.0m	C8	No need		
			U410S-23MR	Adult soft/3.0m	D1			
			U810-23MR	Multi-site Y/3.0m	C8			
DRE								
DRE Waveline EZ Touch Pro			U410-213	Adult finger /3.0m	C7	No need		
			U410S-213	Adult soft /3.0m	C8			
			U810-213	Multi-site Y/3.0m	C7			
Waveline EZ MAX, Waveline Pro (*Oximax Tech)			U410-212P	Adult finger /3.0m	C8	No need		
			U410S-212P	Adult soft/3.0m	D1			
			U810-212P	Multi-site Y/3.0m	C8			
Waveline, Waveline Pro Omni compatible, One			U410-137	Adult finger /3.0m	C7	No need		
			U410S-137	Adult soft/3.0m	C8			
			U810-137	Multi-site Y/3.0m	C7			
Edan (China)								
M9,M99A,M8,M8A, M8B, M3,M3B (*Oximax)		DS100A	U403-01XP	Adult finger /1.1m	C5	U708X-191		
		DS100A	U403S-01XP	Adult soft /1.1m	C7			
		D-YS	U803-01XP	Multi-site Y/1.1m	C5			
EDAN Oximetry M3,M3A,M3B, M8,M8A,M8B M9,M9A,M50,M80 iM8,iM8A,iM8B iM9,iM9A,iM50,iM80 (*Non-Oximax)		MS3-109069	U410-29	Adult finger /3.0m	C7	No need		
			U410S-29	Adult soft /3.0m	C8			
			U810-29	Multi-site Y/3.0m	C7			
		MS-109079	U403-01	Adult finger /1.1m	C3	U708-29		
		DS100A	U403S-01	Adult soft /1.1m	C5			
		D-YS	U803-01	Multi-site Y/1.1m	C3			
M3,M3B,M8,M8A,M8B M9,M9A,M50,M80 iM9,iM9A,iM50,iM80 (*Oximax)			U410-191P	Adult finger /3.0m	C8	No need		
			U410S-191P	Adult soft /3.0m	D1			
			U810-191P	Multi-site Y/3.0m	C8			
(*Digital tech)			U410-46D	Adult finger /3.0m	C8	No need		
			U410S-46D	Adult soft /3.0m	D1			
			U810-46D	Multi-site Y/3.0m	C8			
H100B (*Non Oximax)		MS-109079	U403-01	Adult finger /1.1m	C3	No need		
		DS100A	U403S-01	Adult soft /1.1m	C5			
Edan(China)								

OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.
Edan (China)								
H100N (*Oximax)		DS100A	U403-01XP	Adult finger /1.1m	C5	No need		
		DS100A	U403S-01XP	Adult soft /1.1m	C7			
iT20,iM20, elite V5/V6/V8			U410-48	Adult finger /3.0m	C7	No need		
			U410S-48	Adult soft /3.0m	C8			
			U810-48	Multi-site Y/3.0m	C7			
iT20,iM20, elite V5/V6/V8 (*Nellcor OxiMax)			U410-48P	Adult finger /3.0m	C8	No need		
			U410S-48P	Adult soft /3.0m	D1			
			U810-48P	Multi-site Y/3.0m	C8			
iT20,iM20, elite V5/V6/V8 (Masimo Rainbow)			U410-48MR	Adult finger /3.0m	D1	No need		
			U410S-48MR	Adult soft /3.0m	D2			
			U810-48MR	Multi-site Y/3.0m	D1			
Fukuda Denshi								
DS7100-DS7300 (*Oximax)		DS100A	U403-01XP	Adult finger /1.1m	C5	U710-70		
		DS100A	U403S-01XP	Adult soft /1.1m	C7			
		D-YS	U803-01XP	Multi-site Y/1.1m	C5			
GE Critikon								
Dinamap 9700/9701 9702 PLUS, MPS ORTABLE,/COMPACT, (* Non-Oximax)		DS100A	U403-01	Adult finger /1.1m	C3	U708-01 (EC-8)		
		DS100A	U403S-01	Adult soft /1.1m	C5			
		D-YS	U803-01	Multi-site Y/1.1m	C3			
Dinamap,ProSeries V1 Pro 1000 V2 (Nellcor Optional)		See the above	See the above	See the above	See the above	U708-69 (SCP-10)		
Dinamap,ProCare, ProSeries V2Pro1000 V3 (*Oximax)		DS100A	U403-01XP	Adult finger /1.1m	C5	U710-70 (DOC-10)		
		DS100A	U403S-01XP	Adult soft /1.1m	C7			
		D-YS	U803-01XP	Multi-site Y/1.1m	C5			
Dinamap 8700/8710/8720Pulse, Oxyshuttle11 et 2		DS100A	U403-01	Adult finger /1.1m	C3	Not Available		
		DS100A	U403S-01	Adult soft /1.1m	C5			
		D-YS	U803-01	Multi-site Y/1.1m	C3			
			U410-10	Adult finger /3.0m	C8	No need		
			U410S-10	Adult soft /3.0m	D1			
			U810-10	Multi-site Y/3.0m	C8			
Dinamap ProCare Series V2 Pro 1000 V3 (*M Tech)		1863/ LNCS DCI	U403-49R	Adult finger /1.1m	C4	U708M-15R (2027263-001)		
		1863/ LNCS DCI	U403S-49R	Adult soft /1.1m	C6			
		2258/ LNCS YI	U803-49R	Multi-site Y/1.1m	C4			
		LNOP DC-12	U410-15R	Adult finger /3.0m	C8	No need		
		LNOP DC-12	U410S-15R	Adult soft /3.0m	D1			
			U810-15R	Multi-site Y/3.0m	C8			

OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.
GE Corometrics								
188,129,129F, 128,128F (*Non Oximax)		DS100A	U403-01	Adult finger /1.1m	C3	U708-01		
		DS100A	U403S-01	Adult soft /1.1m	C5			
		D-YS	U803-01	Multi-site Y/1.1m	C3			
120,510,511, 118,556 (*Non Oximax)		See the above	See the above	See the above	See the above	U708-25 (4033CAX)		
128, 129 (*Non Oximax)		See the above	See the above	See the above	See the above	U708-21		
			U410-21	Adult finger /3.0m	C8	No need		
			U410S-21	Adult soft /3.0m	D1			
			U810-21	Multi-site Y/3.0m	C8			
SERIE 250 (*Oximax)		DS100A	U403-01XP	Adult finger /1.1m	C5	U710X-21		
		DS100A	U403S-01XP	Adult soft /1.1m	C7			
		D-YS	U803-01XP	Multi-site Y/1.1m	C5			
		See Oximax Marquette	See Oximax Marquette	See Oximax Marquette	See GE	No need		
GE Datex-Ohmeda								
AS/3,CS/3, AS/5 (Built in M module) CH/S,CH/RSCG/2S, CARDIOCAPI-II Capnomac Ultima Light monitor Satlite, Oscar (Datex tech*)		OXY-F4-N	U410-09	Adult finger /3.0m	C8	No need		
			U110-09	Pediatric finger /3.0m	D1			
			U410S-09	Adult soft /3.0m	D1			
			U110S-09	Pediatric soft /3.0m	D1			
			U210S-09	Infant soft /3.0m	D1			
			U310-09	Neonate wrap /3.0m	D1			
		OXY-W4-N	U610-09	Adult/Neo wrap /3.0m	C8			
			U610B-09	Pedia/infant wrap/3.0m	C8			
		OXY-W4-N	U810-09	Multi-site Y/3.0m	C8			
		OXY-E4-N	U910-09	Adult ear clip /3.0m	D1			
M-NSAT Module (AS/3,CS/3,S/5) (* Non Oximax)		DS100A	U403-01	Adult finger /1.1m	C3	U708-69 (MC-10)		
		DS100A	U403S-01	Adult soft /1.1m	C5			
		D-YS	U803-01	Multi-site Y/1.1m	C3			
M-NSAT Module (i4) (*Oximax)		DS100A	U403-01XP	Adult finger /1.1m	C5	U710-70 (DOC-10)		
		DS100A	U403S-01XP	Adult soft /1.1m	C7			
		D-YS	U803-01XP	Multi-site Y/1.1m	C5			
AS/3,CS/3, AS/5 (Built in M module) CH/S,CH/RSCG/2S, CARDIOCAPI-II Capnomac Ultima Light monitor Satlite, Oscar (Datex tech*)		OXY-F-DB	U403-09	Adult finger /1.1m	C3	U708-09		
			U103-09	Pediatric finger /1.1m	C5			
			U403S-09	Adult soft /1.1m	C5			
			U103S-09	Pediatric soft /1.1m	C5			
			U203S-09	Infant soft /1.1m	C5			
			U303-09	Neonate wrap /1.1m	C5			
		OXY-W-DB	U603-09	Adult/Neo wrap/1.1m	C3			
GE Datex-Ohmeda								


OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.
GE Datex-Ohmeda								
3700,3710,3770, 3775,Biox3740,3760, 3800,3900,3900P, 4700,5250RGB, Cardiacap/5, S/5, AS/3 (built in Module N-XOSAT) CS/3 (Built in module M-OSAT) (*Ohmeda)			U410-02	Adult finger /3.0m	C8	No need		
			U110-02	Pediatric finger /3.0m	D1			
			U410S-02	Adult soft /3.0m	D1			
			U110S-02	Pediatric soft /3.0m	D1			
			U210S-02	Infant soft /3.0m	D1			
			U310-02	Neonate wrap /3.0m	D1			
			U610-02	Adult/Neo wrap /3.0m	C8			
			U610B-02	Pedia/infant wrap/3.0m	C8			
			U810-02	Multi-site Y/3.0m	C8			
			U910-02	Adult ear clip /3.0m	D1			
Ohmeda Tuffsat (*Ohmeda)			U403-02	Adult finger /1.1m	C8	No need		
			U403S-02	Adult soft /1.1m	D1			
Ohmeda TruSat Ohmeda TruSat TD (*+Oxytip+ Technology)			U410-60	Adult finger /3.0m	D1	No need		
			U110-60	Pediatric finger /3.0m	D2			
			U410S-60	Adult soft /3.0m	D2			
			U110S-60	Pediatric soft /3.0m	D2			
			U210S-60	Infant soft /3.0m	D2			
			U310-60	Neonate wrap /3.0m	D2			
			U610-60	Adult/Neo wrap /3.0m	D1			
			U610B-60	Pedia/infant wrap/3.0m	D1			
			U810-60	Multi-site Y/3.0m	D1			
			U910-60	Adult ear clip /3.0m	D2			
			U603B-09	Pedia/infant wrap/1.1m	C3	U708-09		
			U803-09	Multi-site Y/1.1m	C3			
		OXY-E-DB	U903-09	Adult ear clip /1.1m	C5			
GE Hellige								
Cardioserv Eagle1000-3000 (*Non Oximax)		DS100A	U403-01	Adult finger /1.1m	C3	U708-55 (30344358)		
		DS100A	U403S-01	Adult soft /1.1m	C5			
		D-YS	U803-01	Multi-site Y/1.1m	C3			
Cardioserv (* Oximax)		DS100A	U403-01XP	Adult finger /1.1m	C5	NA (2026751-001)		
		DS100A	U403S-01XP	Adult soft /1.1m	C7			
		D-YS	U803-01XP	Multi-site Y/1.1m	C5			
SMU,SMK,SMV (*Ohmeda)		See above *Ohmeda	See above *Ohmeda	See above *Ohmeda	See above *Ohmeda	No need		
GE Marquette								
Dash3000/4000/5000 Solar8000/8000M Solar 8000i/9500 B40i/ 60i/B650 (M Tech)			U410-21MR	Adult finger /3.0m	D2	No need		
			U410S-21MR	Adult soft /3.0m	D3			
			U810-21MR	Multi-site Y/3.0m	D2			
7000/7010,Tram 200/400/600/800 (*Ohmeda Tech.)		See above *Ohmeda	See above *Ohmeda	See above *Ohmeda	See above *Ohmeda			
GE Marquette								

OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext -Cable	Cable Plug	Sensor Conn.
GE Marquette								
Eagle 1000/3000/4000 Dash 2000/3000/4000 Solar 8000/8000M i/9500(Tram 451/451N /851/851N) Cardiolab/Mac-Lab, Corometrics 120/2120 (Tram 451/451N) Responder3000 (nellcor optional) (*Non Oximax)		DS100A	U403-01	Adult finger /1.1m	C3	U708-21 (2006644-001)		
		DS100A	U403S-01	Adult soft /1.1m	C5			
		D-YS	U803-01	Multi-site Y/1.1m	C3			
			U410-21	Adult finger /3.0m	C8	No need		
			U110-21	Pediatric finger /3.0m	D1			
			U410S-21	Adult soft /3.0m	D1			
			U110S-21	Pediatric soft /3.0m	D1			
			U210S-21	Infant soft /3.0m	D1			
			U310-21	Neonate wrap /3.0m	D1			
			U610-21	Adult/Neo wrap/3.0m	C8			
			U610B-21	Pedia/infant wrap/3.0m	C8			
			U810-21	Multi-site Y/3.0m	C8			
			U910-21	Adult ear clip /3.0m	D1			
Dash2500/Dash 3000/ 4000/5000 Solar 8000/8000M i/9500 (Tram 451N5/851N5) Corometrics 250 (GE-Oximax*) Ending "X" can be replaced by "P" for fully compatible with all Oximax oximetry		DS100A	U403-01XP	Adult finger /1.1m	C5	U710X-21 (2021406-001)		
		DS100A	U403S-01XP	Adult soft /1.1m	C7			
		D-YS	U803-01XP	Multi-site Y/1.1m	C5			
			U410-21P	Adult finger /3.0m	D1			
			U110-21P	Pediatric finger /3.0m	D2			
			U410S-21P	Adult soft /3.0m	D2			
			U110S-21P	Pediatric soft /3.0m	D2			
			U210S-21P	Infant soft /3.0m	D2			
			U310-21P	Neonate wrap /3.0m	D2			
			U610-21P	Adult/Neo wrap/3.0m	D1			
			U610B-21P	Pedia/infant wrap/3.0m	D1			
			U810-21P	Multi-site Y/3.0m	D1			
			U910-21P	Adult ear clip /3.0m	D2			
Dash 3000, 4000,5000 Solar 8000/8000M i/9500 (Tram 451M,851M) Cardiolab/Mac-Lab (Tram 451M) Corometrics 259 (*M tech)		1863/ LNCS DCI	U403-49	Adult finger /1.1m	C4	U708M-21R (2027263-002)		
		1863/ LNCS DCI	U403S-49	Adult soft /1.1m	C6			
		2258/ LNCS YI	U803-49	Multi-site Y/1.1m	C4			
			U410-21MR	Adult finger /3.0m	C8	No need		
			U110-21MR	Pediatric finger /3.0m	D1			
			U410S-21MR	Adult soft /3.0m	D1			
			U110S-21MR	Pediatric soft /3.0m	D1			
			U210S-21MR	Infant soft /3.0m	D1			
			U310-21MR	Neonate wrap /3.0m	D1			
			U610-21MR	Adult/Neo wrap/3.0m	C8			
			U610B-21MR	Pedia/infant wrap/3.0m	C8			
			U810-21MR	Multi-site Y/3.0m	C8			
			U910-21MR	Adult ear clip /3.0m	D1			

OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.
GE Ohmeda OxyTip+								
GE S/5 (i4 E series module built in)		OXY-F4-GE	U410-121	Adult finger /3.0m	D1	No need		
			U110-121	Pediatric finger /3.0m	D2			
			U410S-121	Adult soft /3.0m	D2			
			U110S-121	Pediatric soft /3.0m	D2			
			U210S-121	Infant soft /3.0m	D2			
			U310-121	Neonate wrap /3.0m	D2			
			U610-121	Adult/Neo wrap/3.0m	D1			
			U610B-121	Pedia/infant wrap/3.0m	D1			
			U810-121	Multi-site Y/3.0m	D1			
			U910-121	Adult ear clip /3.0m	D2			
Ohmeda TruSat 3755,37XX series 3800,3900,3900P 4700, 5250 Cardicap/5, S/5,AS/3 (N-XOSAT or M built in) CS/3 (M-OSAT built in) S/5 light, Satlite Puls II Capnomac UltmaS/5 ,Coro 259 (*TruSignal tech)		OXY-F-UN	U403-127	Adult finger /1.1m	C8	Cables required depending on the monitor models	Depends on the monitor models See the page of Spo2 adapter cable	
			U103-127	Pediatric finger /1.1m	D1			
		OXY-F-UN	U403S-127	Adult soft /1.1m	D1			
			U103S-127	Pediatric soft /1.1m	D1			
			U203S-127	Infant soft /1.1m	D1			
			U303-127	Neonate wrap /1.1m	D1			
		OXY-W-UN	U603-127	Adult/Neo wrap/1.1m	C8			
		OXY-W-UN	U603B-127	Pedia/infant wrap/1.1m	C8			
		OXY-W-UN	U803-127	Multi-site Y/1.1m	C8			
		OXY-E-UN	U903-127	Adult ear clip /1.1m	D1			
GE Trusignal								
		TS-F-D	U403-117	Adult finger /1.1m	C8	Required	Depends on the monitor models See the page of Spo2 adapter cable	
			U103-117	Pediatric finger /1.1m	D1			
		TS-SA-D	U403S-117	Adult soft /1.1m	D1			
			U103S-117	Pediatric soft /1.1m	D1			
			U203S-117	Infant soft /1.1m	D1			
			U303-117	Neonate wrap /1.1m	D1			
		TS-W-D	U603-117	Adult/Neo wrap/1.1m	C8			
		TS-W-D	U603B-117	Pedia/infant wrap/1.1m	C8			
		TS-W-D	U803-117	Multi-site Y/1.1m	C8			
		TS-E-D	U903-117	Adult ear clip /1.1m	D1			
General (China)								
G1B (*Digital Tech)			U403-42D	Adult finger /1.1m	C5	No need		
			U403S-42D	Adult soft /1.1m	C7			
G2A,G3A,G3C,G3D, G3F,G3G,G3H,G3L, G6A,G6A plus, G9L (*Digital Tech)			U410-65D	Adult finger /3.0m	C8	No need		
			U410S-65D	Adult soft /3.0m	D1			
			U810-65D	Multi-site Y/3.0m	C8			
Guoteng (China)								
PM300		DS100A	U403-01	Adult finger /1.1m	C3	No need		
			U403S-01	Adult soft /1.1m	C5			
			U410-171	Adult finger /3.0m	C8	No need		
			U410S-171	Adult soft /3.0m	D1			
GT3000,GT5000, GT6000GT7000, GT9000,GT9002			U410-179	Adult finger /3.0m	C7	No need		
			U410S-179	Adult soft /3.0m	C8			
			U810-179	Multi-site Y/3.0m	C7			

OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.
Goldway(China)								
UT4000A,UT4000B, UT4000F UT4000ACS			U410-37	Adult finger /3.0m	C7	No need		
			U410S-37	Adult soft /3.0m	C8			
			U810-37	Multi-site Y/3.0m	C7			
		See Nellcor Non-Oximax	See Nellcor Non-Oximax	1.1m	See Nellcor	U708-37		
		See Pace tech	See Pace tech	See Pace tech	See Pace tech	No need		
UT4000A,UT4000C, UT4000EUT4000F, UT4000APro UT4000FPro, UT6000A			U410-11	Adult finger /3.0m	C8	No need		
			U410S-11	Adult soft /3.0m	D1			
			U810-11	Multi-site Y/3.0m	C8			
		See Nellcor Non-Oximax	See Nellcor Non-Oximax	1.1m	See Nellcor	U708-11		
		See Pace tech	See Pace tech	See Pace tech	See Pace tech	No need		
3100		Digital sensor 数字探头	U410-61D	Adult finger /3.0m	C8	No need		
			U410S-61D	Adult soft /3.0m	D1			
			U810-67D	Multi-site Y/3.0m	C8			
UT4000B,UT4000F			U410-67D	Adult finger /3.0m	C8	No need		
			U410S-67D	Adult soft /3.0m	D1			
			U810-67D	Multi-site Y/3.0m	C8			
(*Oximax)			U410-38P	Adult finger /3.0m	D1	No need		
			U410S-38P	Adult soft /3.0m	D2			
			U810-38P	Multi-site Y/3.0m	D1			
UT6000A G3 (*Philips tech)			U410-91	Adult finger /3.0m	C8	No need		
			U410S-91	Adult soft /3.0m	D1			
			U810-91	Multi-site Y/3.0m	C8			
FJBF(China)								
NLB7000N/8000B/8000N NLB_V3,NBL_K3 乐夫:LF-6000A/6000B/ 6000C			U410-134D	Adult finger /3.0m	C8	No need		
			U410S-134D	Adult soft /3.0m	D1			
			U810-134D	Multi-site Y/3.0m	C8			
Infinium Medical								
OmniOne,Omni II Omni III Infitron			U410-97D	Adult finger /3.0m	C8	No need		
			U410S-97D	Adult soft /3.0m	D1			
			U810-97D	Multi-site Y/3.0m	C8			
DRE Waveline Pro (*Digital Tech)			U410-137D	Adult finger /3.0m	C8	No need		
			U410S-137D	Adult soft /3.0m	D1			
			U810-137D	Multi-site Y/3.0m	C8			
DRE Waveline Pro+			U410-147	Adult finger /3.0m	C7	No need		
			U410S-147	Adult soft /3.0m	C8			
			U810-147	Multi-site Y/3.0m	C7			

	OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.
Innocare									
☆	VET-Multi-Parameter			U410-204	Adult finger /3.0m	C7	No need		
				U410S-204	Adult soft /3.0m	C8			
				U810-204	Multi-site Y/3.0m	C7			
IVY									
	405A,405D (* Non Oximax)		DS100A	U403-01	Adult finger /1.1m	C3	U710-85		
			DS100A	U403S-01	Adult soft /1.1m	C5			
			D-YS	U803-01	Multi-site Y/1.1m	C3			
	405P,405T (*Non Oximax)		See the above	See the above	See the above	See the above	U708-01		
	Vital-Guard 450CN(M) (*M tech)		1863/ LNCS DCI	U403-49R	Adult finger /1.1m	C4	U708M-15R		
			1863/ LNCS DCI	U403S-49R	Adult soft /1.1m	C6			
			2258/ LNCS YI	U803-49R	Multi-site Y/1.1m	C4			
			LNOP DC-12	U410-15R	Adult finger /3.0m	C8	No need		
			LNOP DC-12	U410S-15R	Adult soft /3.0m	D1			
				U810-15R	Multi-site Y/3.0m	C8			
☆	450C & CM5 (*Nellcr Oximax)		LNOP DC-12	U410-70P	Adult finger /3.0m	D1	No need		
			LNOP DC-12	U410S-70P	Adult soft /3.0m	D2			
				U810-70P	Multi-site Y/3.0m	D1			
Heal Force									
☆	Classic 90,Classic 120, Prince 100F			U410-61D	Adult finger /3.0m	C8	No need		
				U410S-61D	Adult soft /3.0m	D1			
				U810-61D	Multi-site Y/3.0m	C8			
Jeteem									
	BD-6000 BD-9000			U410-174	Adult finger /3.0m	C7	No need		
				U410S-174	Adult soft /3.0m	C8			
				U810-174	Multi-site Y/3.0m	C7			
JNR									
	DS-100			U410-188	Adult finger /3.0m	C8	No need		
				U410S-188	Adult soft /3.0m	D1			
				U810-188	Multi-site Y/3.0m	C8			
Kontron									
	7840,7845,Minimon 7337Micromon, Supermon,Kolormon, Trakmon (*Non Oximax)			U410-13	Adult finger /3.0m	C8	No need		
				U410S-13	Adult soft /3.0m	D1			
				U810-13	Multi-site Y/3.0m	C8			
			DS100A	U403-01	Adult finger /1.1m	C3	U708-13		
			DS100A	U403S-01	Adult soft /1.1m	C5			
			D-YS	U803-01	Multi-site Y/1.1m	C3			
	7150,7850, Ivy 2000(M)		See the above	See the above	See the above	See the above	U708-01		
			See the above	See the above	See the above	See the above	U708-53		
Lohmeier									
				U410-64	Adult finger /3.0m	D1	No need		
				U410S-64	Adult soft /3.0m	D2			
				U810-64	Multi-site Y/3.0m	D1			

	OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.	
IBIO R&DT										
☆	RD9012			U410-100	Adult finger /3.0m	C7	No need			
				U410S-100	Adult soft /3.0m	C8				
				U810-100	Multi-site Y/3.0m	C7				
☆	(*non Oximax)		DS100A	U410-01	Adult finger /3.0m	C5	U708-56			
			FlexMAX	U410S-01	Adult soft /3.0m	C7				
			D-YS	U810-01	Multi-site Y/3.0m	C5				
*M Technology										
	Rad-8 Rad-5/Rad-5v SatShare		1269/ LNOP DCI	U403-125R	Adult finger /1.1m	C7	Required			
			1276/ LNOP DCIP	U103-125R	Pediatric finger /1.1m	C8				
			1269/ LNOP DCI	U403S-125R	Adult soft 成 /1.1m	C8				
			1276/ LNOP DCIP	U103S-125R	Pediatric soft /1.1m	C8				
				U203S-125R	Infant soft /1.1m	C8				
				U303-125R	Neonate wrap /1.1m	C8				
				U603-125R	Adult/Neo wrap/1.1m	C7				
				U603B-125R	Pedia/infant wrap/1.1m	C7				
				1504/ LNOP YI	U803-125R	Multi-site Y/1.1m				C7
				LNOP TC-I	U903-125R	Adult ear clip /1.1m				C8
			1863/ LNCS DCI	U403-49R	Adult finger /1.1m	C4	Required			
			1864/ LNCS DCIP	U103-49R	Pediatric finger /1.1m	C6				
			1863/ LNCS DCI	U403S-49R	Adult soft /1.1m	C6				
				U103S-49R	Pediatric soft /1.1m	C6				
				U203S-49R	Infant soft /1.1m	C6				
				U303-49R	Neonate wrap /1.1m	C6				
				U603-49R	Adult/Neo wrap/1.1m	C4				
				U603B-49R	Pedia/infant wrap/1.1m	C4				
				2258/ LNCS YI	U803-49R	Multi-site Y/1.1m				C4
				LNCS TC-I	U903-49R	Adult ear clip /1.1m				C6
	M Tech *Refurbished Only		LNOP DCSC	U402-15R	Adult finger /0.9m	C8	No need	Handheld		
			LNOP DC-12	U410-15R	Adult finger /3.0m	C8				
				U110-15R	Pediatric finger /3.0m	D1				
			LNOP DC-12	U410S-15R	Adult soft /3.0m	D1				
				U110S-15R	Pediatric soft /3.0m	D1				
				U210S-15R	Infant soft /3.0m	D1				
				U310-15R	Neonate wrap /3.0m	D1				
				U610-15R	Adult/Neo wrap/3.0m	C8				
				U610B-15R	Pedia/infant wrap/3.0m	C8				
				U810-15R	Multi-site Y/3.0m	C8				
			U910-15R	Adult ear clip /3.0m	D1					
☆	Radical-7,Rad-87 Rad-57,Pronto-7, Rad-8Masimo Rainbow		DCI-dc3	U410-87	Adult finger/3.0m	D2	No need			
				U410S	Adult soft/3.0m	D3				
				U810-87	Multi-site Y/3.0m	D2				

SpO2

OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.
MDE/Invivo								
Escort 1 Escort 300 (*Non Oximax)		DS100A	U403-01	Adult finger /1.1m	C3	U710-85		
		DS100A	U403S-01	Adult soft /1.1m	C5			
		D-YS	U803-01	Multi-site Y/1.1m	C3			
20300,20401,20403, Prism (*Non Oximax)		See the above	See the above	See the above	See the above	U708-01		
Millenia, Escort,(M) M10TM,M12TM, M8TM (*M Tech)		1863/ LNCS DCI	U403-49R	Adult finger /1.1m	C4	U708M-15R		
		1863/ LNCS DCI	U403S-49R	Adult soft /1.1m	C6			
		2258/ LNCS YI	U803-49R	Multi-site Y/1.1m	C4			
		LNOP DC-12	U410-15R	Adult finger /3.0m	C8	No need		
		LNOP DC-12	U410S-15R	Adult soft /3.0m	D1			
4500, ,Omni Track Omeg 1445Millennia			U410-14	Adult finger /3.0m	C8	No need		
			U410S-14	Adult soft /3.0m	D1			
			U810-14	Multi-site Y/3.0m	C8			
4500, ,Omni Track Omeg 1445Millennia M Tech			U410-14M	Adult finger /3.0m	C8	No need		
			U410S-14M	Adult soft /3.0m	D1			
			U810-14M	Multi-site Y/3.0m	C8			
Meditronic - Physio Control								
Lifepack 12 (*non Oximax)		DS100A	U403-01	Adult finger /1.1m	C3	U708-01		
		DS100A	U403S-01	Adult soft /1.1m	C5			
		D-YS	U803-01	Multi-site Y/1.1m	C3			
Lifepack 12, 20(M) (*M Tech)		1863/ LNCS DCI	U403-49R	Adult finger /1.1m	C4	U708M-15R		
		1863/ LNCS DCI	U403S-49R	Adult soft /1.1m	C6			
		2258/ LNCS YI	U803-49R	Multi-site Y/1.1m	C4			
		LNOP DC-12	U410-15R	Adult finger /3.0m	C8	No need		
		LNOP DC-12	U410S-15R	Adult soft /3.0m	D1			
			U810-15R	Multi-site Y/3.0m	C8			
Medlab								
NANOX 10 POX 10 Perarl 10 (*BCI)		3044	U403-06	Adult finger /1.1m	C3	U708-76		
			U403S-06	Adult soft /1.1m	C5			
		3043	U803-06	Multi-site Y/1.1m	C3			
MEK								
MP 100,110, 400,1000			U410-26	Adult finger /3.0m	C6	No need		
			U410S-26	Adult soft /3.0m	C7			
			U810-26	Multi-site Y/3.0m	C6			
(*non Oximax)		DS100A	U403-01	Adult finger /1.1m	C3	U708-56		
		DS100A	U403S-01	Adult soft /1.1m	C5			
		D-YS	U803-01	Multi-site Y/1.1m	C3			
MEK								
			U410-56	Adult finger /3.0m	C8	No need		
			U410S-56	Adult soft /3.0m	D1			
			U810-56	Multi-site Y/3.0m	C8			

OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.			
Mennen											
Envoy (*non Oximax)		DS100A	U403-01	Adult finger /1.1m	C3	U708-54					
		DS100A	U403S-01	Adult soft /1.1m	C5						
		D-YS	U803-01	Multi-site Y/1.1m	C3						
Horizon 1100,2000XL Mercury,MR1300 (*non Oximax)		See the above	See the above	See the above	See the above	U710-85					
Horizon XL,Mercury 1200,MR1300 (*non Oximax)		See the above	See the above	See the above	See the above	U708-24					
Horizon XL (*Ohmeda)		OXY-F4-H	U410-02	Adult finger /3.0m	C8	No need					
			U410S-02	Adult soft /3.0m	D1						
		OXY-W4-H	U810-02	Multi-site Y/3.0m	C8						
Horizon,Envoy, Vitalogik (M) (*M Tech)		1863/ LNCS DCI	U403-49	Adult finger /1.1m	C4	U708M-59					
			U403S-49	Adult soft /1.1m	C6						
		2258/ LNCS YI	U803-49	Multi-site Y/1.1m	C6						
Mindray (China)											
509B, MEC-1000, MEC-1200 New PM6000 PM7000/8000/9000, VS 800,PM600 PM9000 Express (*Mindray Tech)			U410-29	Adult finger /3.0m	C7	No need					
			U110-29	Pediatric finger /3.0m	C8						
			U410S-29	Adult soft /3.0m	C8						
			U110S-29	Pediatric soft /3.0m	C8						
			U210S-29	Infant soft /3.0m	C8						
			U310-29	Neonate wrap /3.0m	C8						
			U610-29	Adult/Neo wrap/3.0m	C7						
			U610B-29	Pedia/infant wrap/3.0m	C7						
			U810-29	Multi-site Y/3.0m	C7						
			U910-29	Adult ear clip/3.0m	C8						
		512D	U403-01	Adult finger /1.1m	C3	U708-29 (0010-20-42594)					
		DS100A	U403S-01	Adult soft /1.1m	C5						
		D-YS	U803-01	Multi-site Y/1.1m	C3						
PM-5000,PM 6000 PM-8000 (*Mindray Tech)		See the above	See the above	See the above	See the above	U708-22 (0515-30-11221)					
			U410-22	Adult finger /3.0m	C8	No need					
			U410S-22	Adult soft /3.0m	D1						
			U810-22	Multi-site Y/3.0m	C8	No need					
			U410-48MR	Adult finger /3.0m	D1				No need		
			U410S-48MR	Adult soft /3.0m	D2						
Beneview T5/T8 (Masimo Rainbow)			U810-48MR	Multi-site Y/3.0m	D1	No need					
Mindray (China)											

Mindray (China)								
OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.
Mindray (China)								
PM-6000 (*M tech)		1863/ LNCS DCI	*U403-49	Adult finger /1.1m	C4	U708M-22 (0010-30-12452)		
		1863/ LNCS DCI	*U403S-49	Adult soft /1.1m	C6			
		2258/ LNCS YI	*U803-49	Multi-site Y/1.1m	C4			
			*U410-22M	Adult finger /3.0m	C8	No need		
			*U410S-22M	Adult soft /3.0m	D1			
New PM-6000, PM-8000 (*M tech)			*U810-22M	Multi-site Y/3.0m	C8	No need		
			*U410-29M	Adult finger /3.0m	C7			
			*U410S-29M	Adult soft /3.0m	C8			
New PM-6000 PM7000,8000,9000 PM-9000 Express VS 800 New PM6000 (*M tech)		See the above M sensor	See the above M sensor	See the above M sensor	See the above	U708M-30		
			*U410-30M	Adult finger /3.0m	C7	No need		
			*U410S-30M	Adult soft /3.0m	C8			
New PM-6000 PM7000,8000,9000 VS 800 PM-9000 Express (*Oximax)		DS100A	U403-01XP	Adult finger /1.1m	C5	U710X-180		
		DS100A	U403S-01XP	Adult soft /1.1m	C7			
		D-YS	U803-01XP	Multi-site Y/1.1m	C5			
			U410-180P	Adult finger /3.0m	C8	No need		
			U410S-180P	Adult soft /3.0m	D1			
Beneview T5/T8 (*Mindray Tech)		See the above Non Oximax Sensor	See the above Non Oximax Sensor	1.1m	See the above	U708-48R		
			U410-48	Adult finger /3.0m	C8	No need		
			U410S-48	Adult soft /3.0m	D1			
			U810-48	Multi-site Y/3.0m	C8			
Beneview T5/T8 (Oximax*)		See Nellcor Oximax	See Nellcor Oximax	1.1m	See Oximax	U710X-48		
Beneview T5/T8 (*M tech)		See M	See M	1.1m	See M	U708M-48		
503,507,508, 509,PM9000 (*Criticare tech)			U410-05	Adult finger /3.0m	C8	No need		
			U410S-05	Adult soft /3.0m	D1			
			U810-05	Multi-site Y/3.0m	C8			
500A,509B, PM9000 (*BCI)			U410-31	Adult finger /3.0m	C8	No need		
			U410S-31	Adult soft /3.0m	D1			
			U810-31	Multi-site Y/3.0m	C8			
PM-50,PM-60 (*Non-Oximax)		DS100A	U403-01	Adult finger /1.1m	C3	No need		
		DS100A	U403S-01	Adult soft /1.1m	C5			

OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.
Minolta								
PULSOX-300		SR-5C	U401-105	Adult finger /0.3m	C8	No need		
PULSOX-300i		SR-5C	U401S-105	Adult soft /0.3m	D1			
M&B(China)								
MB526,MB526+,MB526T8			U410-92	Adult finger /3.0m	C7	No need		
MB526T12,MB526+tw			U410S-92	Adult soft /3.0m	C8			
MB800,CD2000,CD2000+			U810-92	Multi-site Y/3.0m	C7			
CD2000,CD2000+ (* Oximax tech)			U410-92P	Adult finger /3.0m	C8	No need		
			U410S-92P	Adult soft /3.0m	D1			
			U810-92P	Multi-site Y/3.0m	C8			
Nellcor-Tyco								
NPB 40, NPB20, NPB75 NPB 180, NPB 185, NPB-190, NPB-195, NPB 3900, NPB 4000 (*Non-Oximax)+		DS100A	U403-01	Adult finger /1.1m	C3	U704-01 U708-01		
		OXI-P/I	U103-01	Pediatric finger /1.1m	C5			
		DS100A	U403S-01	Adult soft /1.1m	C5			
		OXI-P/I	U103S-01	Pediatric soft /1.1m	C5			
		OXI-P/I	U203S-01	Infant soft /1.1m	C5			
		OXI-A/N	U303-01	Neonate wrap /1.1m	C5			
		D-YS	U803-01	Multi-site Y/1.1m	C3			
		OXI-A/N	U603-01	Adult/Neowrap/1.1m	C3			
		OXI-P/I	U603B-01	Pedia/infantwrap/1.1m	C3			
		D-YS+YSE	U903-01	Adult ear clip /1.1m	C5			
			U410-01	Adult finger /3.0m	C5	No need		
			U110-01	Pediatric finger /3.0m	C7			
			U410S-01	Adult soft /3.0m	C7			
			U110S-01	Pediatric soft /3.0m	C7			
			U210S-01	Infant soft /3.0m	C7			
			U310-01	Neonate wrap /3.0m	C7			
			U610-01	Adult/Neowrap/3.0m	C5			
			U610B-01	Pedia/infantwrap/3.0m	C5			
			U810-01	Multi-site Y/3.0m	C5			
			U910-01	Adult ear clip /3.0m	C7			
N-100,N-200		See the above	See the above	1.1m	See the above	U710-85		
N550,N560,N595, N600,N751, N5500		See Nellcor Oximax	See Nellcor Oximax	1.1m	See Oximax	U710-70		
NPB3000-3100, Symphony, NPB-290, NPB-295, NPB-395		See the above	See the above	1.1m	See the above	U708-69		
Nellcor-Tyco								

↳ Nellcor-Tyco

OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.
Nellcor-Tyco								
			U410-69	Adult finger /3.0m	C8	No need		
			U110-69	Pediatric finger /3.0m	D1			
			U410S-69	Adult soft /3.0m	D1			
			U110S-69	Pediatric soft /3.0m	D1			
			U210S-69	Infant soft /3.0m	D1			
			U310-69	Neonate wrap /3.0m	D1			
			U610-69	Adult/Neowrap/3.0m	C8			
			U610B-69	Pedia/infant wrap/3.0m	C8			
			U810-69	Multi-site Y/3.0m	C8			
			U910-69	Adult ear clip /3.0m	D1			
N-65 (*Oximax) Ending “X” can be replaced by “XP” or “P” for fully compatible with all Oximax oximetry		DS100A	U403-01XP	Adult finger /1.1m	C5	No need		
		OXI-P/I	U103-01XP	Pediatric finger /1.1m	C7			
		OXI-A/N	U403S-01XP	Adult soft /1.1m	C7			
		OXI-P/I	U103S-01XP	Pediatric soft /1.1m	C7			
		OXI-P/I	U203S-01XP	Infant soft /1.1m	C7			
		OXI-A/N	U303-01XP	Neonate wrap /1.1m	C7			
		OXI-A/N	U603-01XP	Adult/Neowrap/1.1m	C5			
		OXI-P/I	U603B-01XP	Pedia/infant wrap/1.1m	C5			
		D-YS	U803-01XP	Multi-site Y/1.1m	C5			
		D-YS+YSE	U903-01XP	Adult ear clip /1.1m	C7			
		See Nellcor Non-Oximax	See Nellcor Non-Oximax	1.1m	See Nellcor	U701-01X		
	Neusoft (China)							
NSC series (*Oximax Tech)			U410-78P	Adult finger /3.0m	C8	No need		
			U410S-78P	Adult soft /3.0m	D1			
			U810-78P	Multi-site Y/3.0m	C8			
			U410-107	Adult finger /3.0m	C7	No need		
			U410S-107	Adult soft /3.0m	C8			
			U810-107	Multi-site Y/3.0m	C7			
New tech (China)								
Pulse Oximer NT1A,NT1B,NT1D (Digital Sensor)			U403-42D	Adult finger /1.1m	C5	No need		
			U403S-42D	Adult soft /1.1m	C7			
Patient monitor Digital sensor			U410-65D	Adult finger /3.0m	C8	No need		
			U410S-65D	Adult soft /3.0m	D1			
			U810-65D	Multi-site Y/3.0m	C8			
NT2A,NT2C			U410-165D	Adult finger /3.0m	D1	No need		
			U410S-165D	Adult soft /3.0m	D2			
			U810-165D	Multi-site Y/3.0m	D1			
NT3A,NT3B, NT3E,NT3F, NT3B60			U410-167	Adult finger /3.0m	C8	No need		
			U410S-167	Adult soft 成/3.0m	D1			
			U810-167	Multi-site Y/3.0m	C8			
NT3A,NT3B, NT3E,NT3F,NT3B60 (*Digital Tech)			U410-173D	Adult finger /3.0m	D1	No need		
			U410S-173D	Adult soft /3.0m	D2			
			U810-173D	Multi-site Y/3.0m	D1			

OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.
Nihon Kohden								
BSM-1101/2101 /23101 7100,7200, OLV-1100 ALife scope			U410-16	Adult finger /3.0m	C8	No need		
			U110-16	Pediatric finger /3.0m	D1			
			U410S-16	Adult soft /3.0m	D1			
			U110S-16	Pediatric soft /3.0m	D1			
			U210S-16	Infant soft /3.0m	D1			
			U310-16	Neonate wrap /3.0m	D1			
			U610-16	Adult/Neowrap/3.0m	D1			
			U610B-16	Pedia/Infant wrap/3.0m	C8			
			U810-16	Multi-site Y/3.0m	C8			
			U910-16	Adult ear clip /3.0m	D1			
Lifescope S, BSS-9800 Lifescope M, BSM 9510 Lifescope P BSM 4100 Life scope A, BSM 5100 Lifescope L, BSM 2350 Lifescope I BSM 2301K, 2303K		TL-201T	U405-16	Adult finger /1.5m	C7	U710-36		
			U105-16	Pediatric finger /1.5m	C8			
			U405S-16	Adult soft /1.5m	C8			
			U105S-16	Pediatric soft /1.5m	C8			
			U205S-16	Infant soft /1.5m	C8			
			U305-16	Neonate wrap /1.5m	C8			
			U605-16	Adult/Neowrap/1.5m	C7			
			U605B-15	Pedia/Infant wrap/1.5m	C7			
			U805-16	Multi-site Y/1.5m	C7			
			U905-16	Adult ear clip /1.5m	C8			
		TL-220T	U410-36	Adult finger /3.0m	D1	No need		
			U110-36	Pediatric finger /3.0m	D2			
			U410S-36	Adult soft /3.0m	D2			
			U110S-36	Pediatric soft /3.0m	D2			
OPV-1500		The same as above	The same as above	1.1m	The same as above	U710-106		
SVM-7501		TL-651S	U410-224	Adult finger /3.0m	C8	No need		
			U410S-224	Adult soft /3.0m	D1			
			U810-224	Multi-site Y/3.0m	C8			

OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.
Nonin								
8500 2500 9840 Avant 4000,9600, 9700 Avant 2120 PulseSense Model 7500,7500FO (Nonin*)		8000AA	U403-08	Adult finger /1.1m	C3	No need		
		8000AP	U103-08	Pediatric finger /1.1m	C5			
		8000SL	U403S-08	Adult soft /1.1m	C5			
		8000SM	U103S-08	Pediatric soft /1.1m	C5			
			U203S-08	Infant soft /1.1m	C5			
			U303-08	Neonate wrap /1.1m	C5			
			U603-08	Adult/Neowrap /1.1m	C5			
			U603B-08	Pedia/infant wrap /1.1m	C3			
			U803-08	Multi-site Y /1.1m	C3			
		8000Q2	U903-08	Adult ear clip /1.1m	C5			
Avant 4000 Wearable Digital Pulse Oximetry			U410-08L	Adult finger /3.0m	C5	No need		
			U401-146	Adult finger /0.3m	D1			
			U401S-146	Adult soft /0.3m	D1			
8600		The same as above	The same as above	1.1m	The same as above	NA		
8604		The same as above	The same as above	1.1m	The same as above	U708-08		
			U410-08	Adult finger /3.0m	C8	No need		
Novametrix								
512/513/710/715, MARS			U403-73	Adult finger /1.1m	C5	No need		
			U103-73	Pediatric finger /1.1m	C7			
			U403S-73	Adult soft /1.1m	C7			
			U103S-73	Pediatric soft /1.1m	C7			
			U203S-73	Infant soft /1.1m	C7			
			U303-73	Neonate wrap /1.1m	C7			
			U603-73	Adult/Neowrap /1.1m	C5			
			U603B-73	Pedia/infant wrap /1.1m	C5			
			U803-73	Multi-site Y /1.1m	C5			
			U903-73	Adult ear clip /1.1m	C7			
			U403-103	Adult finger /1.1m	C8	NA		
			U103-103	Pediatric finger /1.1m	D1			
			U403S-103	Adult soft /1.1m	D1			
			U103S-103	Pediatric soft /1.1m	D1			
			U203S-103	Infant soft /1.1m	D1			
			U303-103	Neonate wrap /1.1m	D1			
			U603-103	Adult/Neowrap /1.1m	C8			
			U603B-103	Pedia/infant wrap /1.1m	C8			
			U803-103	Multi-site Y /1.1m	C8			
			U903-103	Adult ear clip /1.1m	D1			
Novametrix								


Novametrix

OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.
Novametrix								
505/510/511 /515/B/C/520 7000/7100,Oxypleth (Novametrix*)			U410-03	Adult finger /3.0m	C8	No need		
			U110-03	Pediatric finger/3.0m	D1			
			U410S-03	Adult soft /3.0m	D1			
			U110S-03	Pediatric soft/3.0m	D1			
			U210S-03	Infant soft /3.0m	D1			
			U310-03	Neonate wrap /3.0m	D1			
			U610-03	Adult/Neo wrap/3.0m	C8			
			U610B-03	Pedia/infant wrap/3.0m	C8			
			U810-03	Multi-site Y/3.0m	C8			
			U910-03	Adult ear clip /3.0m	D1			
PaceTech/Generra								
		See Nellcor Non-Oximax	See Nellcor Non-Oximax	1.1m	See Nellcor	U708-38		
			U410-38	Adult finger /3.0m	C7	No need		
			U410S-38	Adult soft /3.0m	C8			
			U810-38	Multi-site Y/3.0m	C7			
		See Nellcor Non-Oximax	See Nellcor Non-Oximax	1.1m	See Nellcor	U708-39		
			U410-39	Adult finger /3.0m	C8	No need		
			U410S-39	Adult soft /3.0m	D1			
			U810-39	Multi-site Y/3.0m	C8			
Mediaid/Palco								
M300/M305/340, M5305/M5340			U410-18	Adult finger /3.0m	C8	No need		
			U110-18	Pediatric finger /3.0m	D1			
			U410S-18	Adult soft /3.0m	D1			
			U110S-18	Pediatric soft /3.0m	D1			
			U210S-18	Infant soft /3.0m	D1			
			U310-18	Neonate wrap /3.0m	D1			
			U610-18	Adult/Neo wrap/3.0m	C8			
			U610B-18	Pedia/infant wrap/3.0m	C8			
			U810-18	Multi-site Y/3.0m	D1			
M30/40,M100,M130, M900,M960P, M6 Aster,M4 Nova, M1000			U410-148	Adult finger /3.0m	D1.3	No need		
			U110-148	Pediatric finger /3.0m	D1.6			
			U410S-148	Adult soft /3.0m	D1.6			
			U110S-148	Pediatric soft /3.0m	D1.6			
			U210S-148	Infant soft /3.0m	D1.6			
			U310-148	Neonate wrap /3.0m	D1.6			
			U610-148	Adult/Neo wrap/3.0m	D1.3			
			U610B-148	Pedia/infant wrap/3.0m	D1.3			
			U810-148	Multi-site Y/3.0m	D1.6			

OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.
Philips/HP								
78352C,78354C, 78834C,M1020A, M1025B,M1025A, M1350B,M1772A/B, M1732A/B,M2475B (*Non-Oximax)		M1190A	U410-20	Adult finger /3.0m	C8	No need		
			U110-20	Pediatric finger /3.0m	D1			
		M1190A	U410S-20	Adult soft /3.0m	D1			
			U110S-20	Pediatric soft /3.0m	D1			
			U210S-20	Infant soft /3.0m	D1			
			U310-20	Neonate wrap /3.0m	D1			
			U610-20	Adult/Neo wrap/3.0m	C8			
			U610B-20	Pedia/infant wrap3.0m	C8			
			U810-20	Multi-site Y/3.0m	C8			
			U910-20	Adult ear clip /3.0m	D1			
C1,C3,M3,M4, VirDia, MP20,MP30, MP40/50/60/70/90. VM4,VM6,VM8 (Philips*)		See Philips' sensor	See Philips' sensor	See Philips' sensor	See Philips' sensor	U708-41		
		See Nellcor Non-Oximax	See Nellcor Non-Oximax	1.1m	See Nellcor	U708-43		
A1, A3 (*Non-Oximax)		See Nellcor Non-Oximax	See Nellcor Non-Oximax	1.1m/3.0m	See Nellcor	U708-01		
MP20/30/40/50 PM60/70/90 VN4,VM6,VM8 (*Oximax)		DS100A	U403-01P	Adult finger /1.1m	C5	U710X-43		
		DS100A	U403S-01P	Adult soft /1.1m	C7			
		D-YS	U803-01P	Multi-site Y/1.1m	C5			
MP20/30/40/ PM50/60/70/90 VN4,VM6,VM8 (*M tech)		See M	See M	1.1m	See M	U708M-43R		
			U410-91MR	Adult finger /3.0m	D1	No need		
			U410S-91MR	Adult soft /3.0m	D2			
			U810-91MR	Multi-site Y/3.0m	D1			
MP20/30/40/50/60 70/90/VN4,VM6,VM8 (*Oximax)			U410-91P	Adult finger /3.0m	D1	No need		
			U410S-91P	Adult soft /3.0m	D2			
			U810-91P	Multi-site Y/3.0m	D1			
(Non Oximax)			U410-01	Adult finger /3.0m	C5	U710X-43		
			U410S-01	Adult soft /3.0m	C7			
			U810-01	Multi-site Y/3.0m	C5			
(*Oximax)			U410-01XP	Adult finger /3.0m	C7	U708-41		
			U410S-01XP	Adult soft /3.0m	C8			
			U810-01XP	Multi-site Y/3.0m	C7			

OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.
Philips/HP								
78352C,78354C, 78834C,M1020A, M1025B,M1025A, M1350B,M1772A/B, M1732A/B,M2475B (*Non-Oximax)		DS100A	U403-01	Adult finger /1.1m	C3	U708-20		
		DS100A	U403S-01	Adult soft /1.1m	C5			
		D-YS	U803-01	Multi-site Y/1.1m	C3			
		M1196A	U403-91 U410-91	Adult finger /1.1m 3.0m	C8	U708-40		
		M1192A	U103-91	Pediatric finger /1.1m	D1			
		M1191A	U403S-91 U410S-91	Adult soft /1.1m 3.0m	D1			
		M1192A	U103S-91	Pediatric soft /1.1m	D1			
		M1195A	U203S-91	Infant soft /1.1m	D1			
		M1193A	U303-91	Neonate wrap /1.1m	D1			
			U603-91	Adult/Neo wrap/1.1m	C8			
			U603B-91	Pedia/infant wrap/1.1m	C8			
			U803-91	Multi-site Y/1.1m	D1			
		M1194A	U903-91	Adult ear clip /1.1m	D1			
Schiller ODM/Bruker								
Defigard4000, Defigard5000 ARGUS PRO LifeCare Masimo SET			U410-63M	Adult finger /3.0m	D1	No need		
			U410S-63M	Adult soft /3.0m	D2			
			U810-63M	Multi-site Y/3.0m	D1			
ARGUS PRO LifeCare 2 DEFIGARD Touch 7 Masimo Rainbow			U410-63MR	Adult finger /3.0m	D2	No need		
			U410S-63MR	Adult soft /3.0m	D3			
			U810-63MR	Multi-site Y/3.0m	D2			
Spacelabs								
Ultraview (*Masimo Rainbow)			U410-74MR	Adult finger /3.0m	D2	No need		
			U410S-74MR	Adult soft /3.0m	D3			
			U810-74MR	Multi-site Y/3.0m	D2			
IBIO R&DT								
RD9012			U410-100	Adult finger /3.0m	C7	No need		
			U410S-100	Adult soft /3.0m	C8			
			U810-100	Multi-site Y/3.0m	C7			
ARGUS PRO LifeCare 2 DEFIGARD Touch 7 Masimo Rainbow		DS100A	U410-01	Adult finger /3.0m	C5	No need		
		FlexMAX	U410S-01	Adult soft /3.0m	C7			
		D-YS	U810-01	Multi-site Y/3.0m	C5			
☆Triton								
Defigard4000, Defigard5000 ARGUS PRO LifeCare Masimo SET			U410-135	Adult finger /3.0m	C7	No need		
			U410S-135	Adult soft /3.0m	C8			
			U810-135	Multi-site Y/3.0m	C7			
ARGUS PRO LifeCare 2 DEFIGARD Touch 7 Masimo Rainbow			U410-236	Adult finger /3.0m	C8	No need		
			U410S-236	Adult soft /3.0m	D1			
			U810-236	Multi-site Y/3.0m	C8			

OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.
☆UTAS								
UM-300			U410-129	Adult finger /3.0m	C8			
			U410S-129	Adult soft /3.0m	D1			
			U810-129	Adult ear clip /3.0m	C8			
Unicare(China)								
F8/S			U410-120	Adult finger /3.0m	C7			
			U410S-120	Adult soft /3.0m	C8			
			U810-120	Adult ear clip /3.0m	C7			
Risingmed								
RPM-9000A			U410-217	Adult finger /3.0m	C7	No need		
			U410S-217	Adult soft/3.0m	C8			
			U810-217	Multi-site Y/3.0m	C7			
		DS100A	U403-01	Adult finger /1.1m	C3	U708-217		
		DS100A	U403S-01	Adult soft/1.1m	C5			
		D-YS	U803-01	Multi-site Y/1.1m	C3			
RGB medical								
(*BCI)		3044	U403-06	Adult finger /1.1m	C3	U708-76		
			U403S-06	Adult soft /1.1m	C3			
		3078	U803-06	Adult ear clip /1.1m	C3			



OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext -Cable	Cable Plug	Sensor Conn.
Schiller ODM/Bruker								
Propaq CS, Encore Series, Atlas (* non-Oximax)		DS100A	U403-01	Adult finger /1.1m	C3	U708-01		
		DS100A	U403S-01	Adult soft /1.1m	C5			
		D-YS	U803-01	Multi-site Y/1.1m	C3			
TM 910-Defigard 500(M) (* non-Oximax)		See Nellcor Non-Oximax	See Nellcor Non-Oximax	1.1m/3.0m	See Nellcor	U708-63		
Argus VCM (*Oximax)		DS100A	U403-01XP	Adult finger /1.1m	C5	U710-70		
		DS100A	U403S-01XP	Adult soft /1.1m	C7			
		D-YS	U803-01XP	Multi-site Y/1.1m	C5			
CM8 SM783,784,785, TM910		See RGB medical	See RGB medical	See RGB medical	See RGB medical	U708-17		
			U410-17	Adult finger /3.0m	C8	No need		
			U410S-17	Adult soft /3.0m	D1			
TM7,Argus LCM, (*M tech)		1863/ LNCS DCI	U403-49	Adult finger /1.1m	C4	U708M-17		
		1863/ LNCS DCI	U403S-49	Adult soft /1.1m	C6			
		2258/ LNCS YI	U803-49	Multi-site Y/1.1m	C4			
			U410-17M	Adult finger /3.0m	C8	No need		
			U410S-17M	Adult soft /3.0m	D1			
			U810-17M	Multi-site Y/3.0m	C8			
Propaq CS, Encore, Micropaq (New) (M) (*M tech)		See M	See M	1.1m	See M	U708M-15R		
		LNOP DC-12	U410-15R	Adult finger /3.0m	C8	No need		
		LNOP DC-12	U410S-15R	Adult soft /3.0m	D1			
				U810-15R	Multi-site Y/3.0m			

OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.
Spacelabs								
90419/01,90651A/07, 90501, Ultraview 90496 (*Novamatrix)			U410-03	Adult finger /3.0m	C8	No need		
			U110-03	Pediatric finger /3.0m	D1			
			U410S-03	Adult soft /3.0m	D1			
			U110S-03	Pediatric soft /3.0m	D1			
			U210S-03	Infant soft /3.0m	D1			
			U310-03	Neonate wrap /3.0m	D1			
			U610-03	Adult/Neo wrap/3.0m	C8			
			U610B-03	Pedia/infant wrap/3.0m	C8			
			U810-03	Multi-site Y/3.0m	C8			
			U910-03	Adult ear clip /3.0m	D1			
90351-0/6,90465, 90466, 90456,PC exoress (*Non-Oximax)		DS100A	U403-01	Adult finger /1.1m	C3	U708-04		
		DS100A	U403S-01	Adult soft /1.1m	C5			
		D-YS	U803-01	Multi-site Y/1.1m	C3			
			U410-04	Adult finger /3.0m	C8	No need		
			U110-04	Pediatric finger /3.0m	D1			
			U410S-04	Adult soft /3.0m	D1			
			U110S-04	Pediatric soft /3.0m	D1			
			U210S-04	Infant soft /3.0m	D1			
			U310-04	Neonate wrap /3.0m	D1			
			U610-04	Adult/Neo wrap/3.0m	C8			
			U610B-04	Pedia/infant wrap/3.0m	C8			
			U810-04	Multi-site Y/3.0m	C8			
			U910-04	Adult ear clip /3.0m	D1			
Ultraview (*Non-Oximax)		See Nellcor Non-Oximax	See Nellcor Non-Oximax	1.1m	See Nellcor	U708-74		
			U410-74	Adult finger /3.0m	D1	No need		
			U410S-74	Adult soft /3.0m	D2			
Ultraview (M Tech)			U810-74	Multi-site Y/3.0m	D1	No need		
			U410-74MR	Adult finger /3.0m	D1			
			U410S-74MR	Adult soft /3.0m	D2			
Ultraview (Oximax Tech)			U810-74MR	Multi-site Y/3.0m	D1	No need		
			U410-74P	Adult finger /3.0m	D1			
			U410S-74P	Adult soft /3.0m	D2			
PRO2 (*BCI)			U810-74P	Multi-site Y/3.0m	D1	No need		
		3044	U403-06	Adult finger /1.1m	C3			
			U403S-06	Adult soft /1.1m	C5			
		3043	U803-06	Multi-site Y/1.1m	C3			

OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.
Takaoka								
Takaoka 9504, 9504ECG, Oxipani , Oxifast			U410-99	Adult finger /3.0m	C7	No need		
			U410S-99	Adult soft /3.0m	C8			
			U810-99	Multi-site Y/3.0m	C7			
Uray								
UR-8000 UR-9000			U410-143	Adult finger /3.0m	C7	No need		
			U410S-143	Adult soft /3.0m	C8			
			U810-143	Multi-site Y/3.0m	C7			
Zoll								
M series AED (*M tech)		1863/ LNCS DCI	U403-49R	Adult finger /1.1m	C4	U708M-15R		
		1863/ LNCS DCI	U403S-49R	Adult soft /1.1m	C6			
		2258/ LNCS YI	U803-49R	Multi-site Y/1.1m	C4			
		LNOP DC-12	U410-15R	Adult finger /3.0m	C8	No need		
		LNOP DC-12	U410S-15R	Adult soft /3.0m	D1			
			U810-15R	Multi-site Y/3.0m	C8			
Zondan (China)								
Apollo N2, Apollo N4 Apollo N5 (M tech optional) DS120C			U410-45	Adult finger /3.0m	C7	No need		
			U410S-45	Adult soft /3.0m	C8			
			U810-45	Multi-site Y/3.0m	C7			
ZD120B ZD120D ZD120E (*Digital Tech)			U410-66D	Adult finger /3.0m	C8	No need		
			U410S-66D	Adult soft /3.0m	D1			
			U810-66D	Multi-site Y/3.0m	C8			
ZD120D/ZD120E			U410-45	Adult finger /3.0m	C7	No need		
			U410S-45	Adult soft /3.0m	C8			
			U810-945	Multi-site Y/3.0m	C7			
Apollo N1/N2/N3/N4/N5			U410-210M	Adult finger /3.0m	D1	No need		
			U410S-210M	Adult soft /3.0m	D2			
			U810-210M	Multi-site Y/3.0m	D1			
Unicare(China)								
(*Digital tech)			U410-66D	Adult finger /3.0m	C8	No need		
			U410S-66D	Adult soft /3.0m	D1			
			U810-66D	Multi-site Y/3.0m	C8			

Please be Reminded : The other types of sensors like U100, U100S, U200S, U300, U600, U600B and U900 of the listed respected brands are available upon demands.