














③ Patient NIBP Adapter Air Hose




Original Maker & Models	Patient Size	Plug to Monitor	Plug to Cuff	Hose Length	REF#	Original REF#	Package	Unit Price
Philips All models except A1/A3	Adult Pediatric			1.5m 2.5m 3.0m	HS15-16-15 HS-16-15 HS30-16-15	M1599B M1598B	1 pcs	B7
	Infant Neonate			1.5m 2.5m 3.0m	HSN15-16-20 HSN-16-20 HSN30-16-20	M1596B/C M1597B/C	1 pcs	C1
Philips A1/A3 Nellcor, Siemens, Datascope, Spacelabs, Kontron, Mindray	Adult Pediatric			2.5m 3.0m	HS-15-15 HS30-15-15	M3918A(Philips)	1 pcs	B7
	Infant Neonate			2.5m	HSN-15-20		1 pcs	C1
GE-Dinamap Pro, Compact, MPS Pro Care series	Adult Pediatric			2.5m 3.6m 7.2m	HD-25-09 HD36-25-09 HD72-25-09	107365 107366	1 pcs	C4
	Infant Neonate			2.5m 3.6m	HDN-25-27 HDN36-25-27	107368	1 pcs	C4
	Adult Pediatric			2.5m 3.6m	HD-25-17 HD36-25-17	107363	1 pcs	C3
GE-Datex-Ohmeda AS/3, CS/3, S/5 Cardiacap 5/II, Light	Adult Pediatric			1.8m 2.5m 3.0m 6.0m	HD18-22-17/18 HD-22-17/18 HD30-22-17/18 HD60-22-17/18	895732 877235 879739	1 pcs	C3
	Infant Neonate			2.5m 3.0m	HDN-22-27 HDN30-22-27	9461-205 877514	1 pcs	C4
	Infant Neonate			3.6m	HDN36-22-20	9461-215	1 pcs	C3
GE-Marquette Eagle4000, Solar Tram100/200/300	Adult Pediatric			2.5 3.6m	HD-22-17 HD36-22-17	9461-203	1 pcs	C3
	Adult Pediatric			2.5m 3.6m	HD-22-09 HD36-22-09	9461-214	1 pcs	C3
	Infant Neonate			2.5m 3.0m	HDN-22-27 HDN30-22-27	9461-205 877514	1 pcs	C4
GE-Marquette -Dash, Eagle, SOLAR Tram400/450/451 800/850/851, Coro 120	Adult Pediatric			2.5m 3.6m	HD-24-17 HD36-24-17	414873-001	1 pcs	C3
	Adult Pediatric			2.5m 3.6m	HD-24-09 HD36-24-09	9461-217	1 pcs	C3
	Infant Neonate			2.5m	HDN-24-27	414874-001	1 pcs	C4
	Infant Neonate			3.6m	HDN36-24-20	9461-216	1 pcs	C3
	Infant Neonate			2.5m	HDN-24-20-2	9461-218	1 pcs	C3


Original Maker & Models	Patient Size	Plug to Monitor	Plug to Cuff	Hose Length	REF#	Original REF#	Package	Unit Price
GE-Critikon Invivo, Kontron, Nihon Kohden	Adult Pediatic			2.5m 3.6m 7.2m	HD-09-09 HD36-09-09 HD72-09-09	008841 08842	1 pcs	C3
	Adult Pediatic			3.6m	HD36-09-17	088845	1 pcs	C3
	Infant Neonate			2.5m	HDN-09-27	008840	1 pcs	C4
Nihon Kohden BSM/PVB	Adult Pediatic			1.5m 2.5m 3.5m	HD15-30-18 HD-30-18 HD35-30-18	YN-900P YN-901P	1 pcs	C6
	Infant Neonate			1.5m 2.5m 3.5m	HDN15-30-20 HDN-30-20 HDN35-30-20	YN-920P YN-921P	1 pcs	C7
	Life Scope N OPV-1500 Life Scope N OPV-1500K			1.5m 2.5m 3.0m 3.5m	HD15-37-18 HD-37-18 HD30-37-18 HD35-37-18	YN-100T YN-101T YN-101T	1 pcs	C6
Criticare, Datascope, Spacelabs, Colin, Welch Allyn, Nihon Kohden	Adult Pediatic			2.5m	HS-04-07		1 pcs	B7
	Infant Neonate			2.5m	HSN-04-20		1 pcs	C1
Colin	Adult Pediatic			2.5m 3.5m	HS-31-32 HS35-31-32	A015ZZ	1 pcs	C1
	Infant Neonate			2.5m 3.5m	HSN-31-20 HSN35-32-20	AL012Z	1 pcs	C1
Drager/Siemens SC6000- 7000/8000/ 9000 Infinity	Adult Pediatic			2.5m 3.7m	HS-40-15 HS37-40-15	MP00953 (1275 275 EH50u)	1 pcs	C6
	Infant Neonate			2.4m 2.5m	HSN24-40-20 HSN-40-20	2870298	1 pcs	C7
Physio Control Lifepack 20	Adult Pediatic			2.5m	HS-12-15		1 pcs	B7
	Infant Neonate			2.5m	HSN-12-20		1 pcs	C1
MEK MP700, MP1000	Adult Pediatic			2.5m	HS-09-09		1 pcs	B7
	Infant Neonate			2.5m	HSN-09-20		1 pcs	C1

Original Maker & Models	Patient Size	Plug to Monitor	Plug to Cuff	Hose Length	REF#	Original REF#	Package	Unit Price
Medtronic	Infant Neonate				HTS-44-20		1 pcs	C8
	Adult Pediatric			2.5m	HS-04-09		1 pcs	B7
Welch Allyn Atlas	Infant Neonate			2.5m	HSN-04-20		1 pcs	C1
	Adult Pediatric			2.5m	HD-41-09		1 pcs	C3
	Infant Neonate			2.5m	HDN-41-27		1 pcs	C3
	Adult Pediatric			2.5m	HS-42-43		1 pcs	C3
Infinium OMNI/OMNI EXPRESS/ OMNI EXPRESS VS OMNI II/OMNI III	Adult Pediatric			2.1m	HS-36-2T		1 pcs	B7






④ NIBP Connectors & Tubing




4.1 Universal Y Adapter

4.1a 1-Tube Air Hose To 2-Tube BP Cuff ,10pcs/pack



REF#	BP3T-Y	BP3T-1718	BP3T-1717	BP3T-0909	BP3T-2T2T	BP3T-2020
Image						
Price	B7	C6	C6	D1	C2	D1
REF#	BP12-3T	BP12-1718	BP12-1717	BP12-0909	BP12-2T2T	BP12-2020
Image						
Price	C2	C8	C8	D3	C6	D3
REF#	BP18-3T	BP18-1718	BP18-1717	BP18-0909	BP18-2T2T	BP18-2020
Image						
Price	C2	C8	C8	D3	C6	D3
REF#	BP04-3T	BP04-1718	BP04-1717	BP04-0909	BP04-2T2T	BP04-2020
Image						
Price	C8	D3	D3	D4	D1	D4
REF#	BP08-3T	BP08-1718	BP08-1717	BP08-0909	BP08-2T2T	BP08-2020
Image						
Price	C2	C8	C8	D3	C6	D3
REF#	BP2T-3T	BP2T-1718	BP2T-1717	BP2T-0909	BP2T-2T2T	BP2T-2020
Image						
Price	C2	C8	C8	D3	C6	D3
REF#	BP21-3T	BP21-1718	BP21-1717	BP21-0909	BP21-2T2T	BP21-2020
Image						
Price	C2	C8	C8	D3	C6	D3

4.1b 2-Tube Air Hose To 1-Tube BP Cuff,10 pcs/pack



REF#	BP3T-Y	BP15-3T	BP04-3T	BP18-3T	BP2T-3T	BP20-3T
Image						
Price	B7	D1	C8	C2	C2	C6
REF#	BP3T-1718	BP15-1718	BP04-1718	BP18-1718	BP2T-1718	BP20-1718
Image						
Price	C6	D4	D3	C8	C8	D1
REF#	BP3T-1818	BP15-1818	BP04-1818	BP18-1818	BP2T-1818	BP20-1818
Image						
Price	C6	D4	D3	C8	C8	D1
REF#	BP3T-0808	BP15-0808	BP04-0808	BP18-0808	BP2T-0808	BP20-0808
Image						
Price	C2	D3	D1	C6	C6	C8
REF#	BP3T-1212	BP15-1212	BP04-1212	BP18-1212	BP2T-1212	BP20-1212
Image						
Price	C6	D4	D3	C8	C8	D1
REF#	BP3T-0505	BP15-0505	BP04-0505	BP18-0505	BP2T-0505	BP20-0505
Image						
Price	D1	D4	D3	D3	D3	D3

4.2 Connectors for Cuffs and Hoses



Unit:Per pcs

REF#	BP03	BP04	BP05	BP06	BP07	
Fig.						
Match	(BP04)>> BP06	(BP03)>>BP07	>>BP20	(BP07)>>BP03	(BP06)>>BP04	
Comp	Mennen (Old) Schiller Bruker/Odam Spacelabs(Old)	Same as left	CAS, Drager Kontron(Amphora) Physio control(Old) Schiller Bruker/Odam	BCI Monitor (Old) Datascope(Old)	Same as left	
Price	A1		A5	A3	A1	A5
REF#	BP08	BP09	BP10	BP12	BP12-P	BP15
Fig.						
Match	>> BP09	>>BP08	Match	(BP12-P)>>BP15	(BP12)>>BP15	>>BP12
Comp	Welch Allyn, GE Critikon(2) Invivo(2), Kontron(2) Mennen(2) Nihon Kohden(2)	GE-Critikon Invivo, Kontron, Nihon Kohden		Criticare, Datascope Drager/Siemens Physio control Mindray,,Philips/HP Spacelabs, Nellcor	Same as BP12	Schiller Air hose connector
Price	A1	A3	A6	A2	A1	A6

Unit:Per pcs





REF#	BP16	BP17	BP18	BP20	BP21	BP2T-M
Fig.						
Match		>>BP18	>>BP17	>>BP05	>>BS21	Two end thread
Comp	Philips	Nihon Kohden Datex(2)	Zoll Series M Datex(1) Marquette(2)	Regularly as the connector of air hose for neonate disposable BP cuffs	Pace Tech	General purpose
Price	A5	A2	A2	A3	A2	A2
REF#	BP22	BP24	BP25	BP26	BP27	BP28
Fig.						
Match	BS22					
Comp	GE-Datex-Ohmeda	GE-Marquette	GE-Dinamap			
Price	A3	A5	A5	A5	A3	A3
REF#	BP30	BP30A	BP31	BP32	BP34	BP35
Fig.						
Match			>.BP35	>>BP31	>.BP35	>>BP31
Comp	Nihon Kohden	Nihon Kohden	Colin Gambro Tyco(NPB 5500)	Colin	Colin Gambro Tyco(NPB 5500)	Monitor plug
Price	A4	A4	A5	A6	A5	A5

Unit:Per pcs

REF#	BP36	BP37	BP38	BP39	BP40	BP41
Fig.						
Match	>>BS36		>>BP39	>>BP38	>>BP12	>>BS41
Comp		Nihon Kohden	Kendall	Kendall	Drager/Siemens	Welch Allyn Altas
Price	A7	A6	A3	A1	A7	A3
REF#	BP42	BP43	BP44	BP45	BP46	BP47
Fig.						
Match	>>BP43	>>BP42				
Comp	Welch Allyn	Welch Allyn	Medtronic	Welch Allyn	Welch Allyn	GE
Price	A4	A3	A6	A2	A2	A4
REF#	BP48	BP49	BP51	BP52	BP54	
Fig.						
Match	>>BP47	>>BP48	>>BP54	>>BP51	>>BP51	
Comp	GE	GE	GE Nihon Kohden SVM-7500	GE Nihon Kohden SVM-7500	GE Nihon Kohden SVM-7500	
Price	A2	A4	A1	A1	A3	

4.3 Air Tubing



Figures	REF#	Application	Inner Diameter	Color	Price / m
	Y02-H008	Adult/Pediatric single tube	Ø3.50mm	Light gray	A3
	Y02-H007	Adult/Pediatric double tube	Ø3.50mm	Light gray	A4
	Y02-H010	Neonate/ Infant single tube	Ø2.70mm	Light blue	A4
	Y02-H009	Neonate/ Infant double tube	Ø2.70mm	Light blue	A6